



MEL/850

Potable Water Laboratory

CERTIFICATION

Hach Company certifies the Portable Incubator was tested thoroughly, inspected and found to meet its published specifications when it was shipped from the factory. The **Portable Incubator** has been tested and is certified as indicated to the following instrumentation standards:

Product Safety: (Power Supplies Only)

Power Supply, 115 or 220 VAC Input, 12 VDC Output: Unlisted

Power Supply, 230 VAC Input, 12 VDC Output: CE/VDE Approved, GS

**Immunity (Instrument Tested with external 230 VAC Power Supply):
EN 50082-1** (European Generic Immunity Standard) **per 89/336/EEC**

EMC: Supporting test records by Hach Company, certified compliance by Hach Company.

Required Standard/s include:

EN 61000-4-2 "1995" (IEC 1000-4-2) Electro-Static Discharge

EN 61000-4-4 "1995" (IEC 1000-4-4) Electrical Fast Transients/Burst

EN 61000-4-11 "1994" (IEC 1000-4-11) Voltage Dips, Interruptions and Variations

ENV 50140 "1993" (IEC 1000-4-3) Radiated RF Electro-Magnetic Fields

ENV 50141 "1993" Conducted Disturbances Induced by RF Fields

ENV 50204 "1995" Radiated Electro-Magnetic Field from Digital Telephones

EN 61000-4-5 "1995" (IEC 1000-4-5) Surge, Light Industrial Levels

Emissions (Instrument Tested with external 230 VAC Power Supply):

EN 50081-1 (Emissions) **per 89/336/EEC EMC:** Supporting test records by Criterion Tech certified compliance by Hach Company.

Required Standard/s include:

EN 55022 (CISPR 22) Emissions, Class A Limits (Testing to certify to "B" Limits by 31 January 1998)

Additional Standard/s include:

EN 61000-2 Harmonic Disturbances Caused by Electrical Equipment

EN 61000-3 Voltage Fluctuation (Flicker) Disturbances Caused by Electrical Equipment

**CANADIAN INTERFERENCE-CAUSING EQUIPMENT
REGULATION, IECS-003, Class A:**

Supporting test records by Criterion Technology, certified compliance by Hach Company.

This Class A digital apparatus meets all requirements of the Canadian Interference-Causing Equipment Regulations.

Cet appareil numérique de la classe A respecte toutes les exigences du Règlement sur le matériel brouilleur du Canada.

CERTIFICATION, continued

FCC PART 15, Class “A” Limits:

Supporting test records by Criterion Technology, certified compliance by Hach Company.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

(1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference, in which case the user will be required to correct the interference at his own expense. The following techniques of reducing the interference problems are applied easily.

1. Disconnect the Portable Incubator from its power source to determine if the Portable Incubator is the source of the interference.
2. If the Portable Incubator's AC power supply is plugged into the same outlet as the device with which it is interfering, try another outlet.
3. Move the Portable Incubator away from the device receiving the interference.
4. Reposition the receiving antenna for the device receiving the interference.
5. Try combinations of the above.

CERTIFICATION, continued

Hach Company certifies the Pocket Pal™ instruments were tested thoroughly, inspected and found to meet their published specifications when it was shipped from the factory. The instruments have been tested and are certified as indicated to the following instrumentation standards:

EN 55011 per 89/336/EEC EM. Tested by Blackstone Inc., certified by Hanna.

EN 50082-1 (Immunity) per 89/336/EEC. Tested by Blackstone Inc., certified by Hanna. Standards include:

IEC 801-2 (ESD)

IEC 801-3 (RF Radiated)

TABLE OF CONTENTS

CERTIFICATION.....	iii
SAFETY PRECAUTIONS	ix
OPERATION.....	1
SECTION 1 INTRODUCTION	3
1.1 Portable Incubator	3
1.2 DR/850 Colorimeter.....	4
1.3 Capabilities of the MEL/850 Potable Water Laboratory	5
1.4 Preparation for Use.....	5
1.5 Missing or Damaged Product.....	6
SECTION 2 PORTABLE INCUBATOR	7
2.1 Specifications	7
2.2 Power Selection.....	7
2.2.1 Using In A Motorized Vehicle	7
2.2.2 Using Portable Battery Pack.....	7
2.2.3 Using Battery Eliminator.....	8
2.3 Operation of Controls and Indicators.....	9
2.3.1 The Display.....	9
2.4 Operation.....	10
2.5 Calibration.....	11
2.6 Maintenance	12
2.6.1 Cleaning.....	12
2.6.2 Replacing the Fuses	12
2.6.2.1 Replacing the 115/230 Volt Battery Eliminator Fuses	13
2.6.2.2 Replacing the Battery Pack Fuse	13
SECTION 3 MICROBIOLOGICAL PROCEDURES	15
3.1 Overview of Coliform Bacteria.....	15
3.2 Testing Techniques	16
3.3 Preparing Sample Containers.....	16
3.3.1 Presterilized Containers.....	16
3.3.2 Autoclavable Containers.....	16
3.4 Collecting and Preserving Samples.....	17
3.4.1 Samples from Faucets, Spigots, Hydrants, or Pumps.....	18
3.4.2 Samples from Rivers, Lakes, and Reservoirs	18
3.5 Disposing of Completed Tests	18
3.6 Presence/Absence Methods—Bromcresol Purple Broth (Method 8319) and Bromcresol Purple Broth with MUG (Method 8364).....	19

TABLE OF CONTENTS, continued

3.6.1 Preparing Materials	19
3.6.2 Convenient Packaging	19
3.7 Using P/A Broth Disposable Bottles	20
3.8 Using P/A Broth Ampules	21
3.9 Interpreting P/A Results.....	22
3.10 Confirming Positive Samples	22
3.11 Summary of Method 8319	23
3.12 Summary of Method 8364	23
SECTION 4 MEASURING TOTAL DISSOLVED SOLIDS (TDS)	25
4.1 Calibrating the TDS Pocket Pal Tester	25
4.2 Replacing the Batteries	26
4.3 TDS Procedure.....	27
SECTION 5 MEASURING pH	29
5.1 Calibrating the pH Pocket Pal Tester.....	29
5.2 Replacing the Batteries	30
5.3 pH Procedure	31
DR/850 COLORIMETER PROCEDURES	33
CHLORINE, FREE, DPD Method	35
CHLORINE, TOTAL, DPD Method.....	45
NITRATE, Low Range, Cadmium Reduction Method	55
NITRITE, Low Range, Diazotization Method	63
NITROGEN, AMMONIA, Salicylate Method.....	69
pH, Colorimetric pH Determination Using Phenol Red	75
PHOSPHORUS, REACTIVE, PhosVer 3 (Ascorbic Acid) Method.....	79
SULFIDE, Methylene Blue Method.....	87
GENERAL INFORMATION.....	91
HOW TO ORDER.....	93
REPAIR SERVICE.....	94
LIMITED WARRANTY	95

SAFETY PRECAUTIONS

Before attempting to unpack, set up, or operate any instruments in this kit, please read the instruction manuals shipped with them. Pay particular attention to all warnings, cautions and notes. Failure to do so could result in serious injury to the operator or damage to the equipment.

To ensure the protection provided by this equipment is not impaired, this equipment **MUST NOT** be installed or used in any manner other than that which is specified in this manual.

Use of Hazard Information

If multiple hazards exist, the signal word corresponding to the greatest hazard shall be used.

DANGER

Indicates either a potentially or an imminently hazardous situation which, if not avoided, could result in either death or serious injury

CAUTION

Indicates a potentially hazardous situation that may result in minor or moderate injury

NOTE

Information that requires special emphasis

Precautionary Labels

Please pay particular attention to labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed.

The DR/800 Series Colorimeters are Class 1 LED products. A Class 1 LED product has insufficient energy to be considered an eye hazard.



This symbol, if noted on the instrument, references the Instruction Manual for operational and/or safety information.



OPERATION

DANGER

Handling chemical samples, standards, and reagents can be dangerous. Review the necessary Material Safety Data Sheets and become familiar with all safety procedures before handling any chemicals.

PERIGO

A manipulação de amostras, padrões e reagentes químicos pode ser perigosa. Reveja as necessárias Fichas Técnicas de Segurança do Material e familiarizese com os procedimentos de segurança antes de manipular quaisquer substâncias químicas.

PELIGRO

La manipulación de muestras químicas, patrones y reactivos puede ser peligrosa. Antes de manipular cualquier productor químico, conviene leer las Fichas Técnicas de Seguridad y familiarizarse con los procedimientos de seguridad.

DANGER

La manipulation des échantillons chimiques, étalons et réactifs peut être dangereuse. Lire les fiches de données de sécurité des produits nécessaires et se familiariser avec toutes les procédures de sécurité avant de manipuler tout produit chimique.

GEFAHR

Da das Arbeiten mit chemikalischen Proben, Standards, Reagenzien und Abfällen mit Gefahren verbunden ist, empfiehlt die Hach Company dem Benutzer dieser Produkte dringend, sich vor der Arbeit mit sicheren Verfahrensweisen und dem richtigen Gebrauch der Chemikalien oder Biogefahrgut vertraut zu machen und alle entsprechenden Material Sicherheitsdatenblätter aufmerksam zu lesen.

SECTION 1 INTRODUCTION

1.1 Portable Incubator

The Hach Portable Incubator is a bacterial incubator designed for field use in Hach's Microbiological Environmental Laboratories (MELs). Available MELs include the MEL P/A Safe Drinking Water Laboratory, the MEL/MPN Laboratory for Total Coliforms and *E. coli*, the MEL/MF Laboratory for Total Coliforms, and the MEL/850 Potable Water Laboratory.

The Portable Incubator maintains temperatures with ± 0.5 °C and the incubation temperature is adjustable between 30 and 50 °C. Ideally suited for total coliform, fecal coliform, and *E. coli* testing, the incubator may be used for Presence/Absence (P/A), Membrane Filtration (MF), and the Most Probable Number (MPN) procedures.

The instrument power cord easily plugs into an automobile cigarette lighter. For remote field use, a 12 Vdc portable battery is available. The portable battery is rechargeable and includes recharger and nylon carrying case. Battery eliminators are also available for 115 and 230 Vac use.

Optional accessories include:

- Rack for P/A bottles (holds 6)
- Rack for MPN tubes (holds 39)
- Rack for MF petri dishes (holds 42 50-mm dishes)
- Portable battery
- Portable eliminators (115 or 230 Vac)
- Testing media and apparatus

See *SECTION 2 PORTABLE INCUBATOR* for operating instructions.

1.2 DR/850 Colorimeter

The Hach DR/850 Colorimeter is a microprocessor-controlled filter photometer with an LED light source. It is suitable for colorimetric testing in the laboratory or the field. The instrument is precalibrated for common colorimetric measurements and includes convenient calibration capability for user-entered and future Hach methods.

Some of the features it offers are:

- Displays results in concentration, absorbance, or % transmittance.
- Automatic wavelength selection and ranging in the preprogrammed parameters.
- Data storage and recall for datalogging in the field or laboratory.
- Icon prompts displayed during testing.
- IR output for RS232 interface capability allows an external printer or computer to interface with the colorimeter.
- Entry of user-entered methods or new Hach methods.
- Error signals for procedural or instrument troubleshooting.

The instrument holds four AA-size alkaline dry cells (batteries supplied) that power the instrument for at least six months. Optional rechargeable alkaline batteries are also available. The charger and optional rechargeable batteries must be purchased separately.

See the *DR/800 Series Instrument Manual* for more information on features, operation, calibration, and maintenance.

SECTION 1, continued

1.3 Capabilities of the MEL/850 Potable Water Laboratory

Parameter	Reasons for Measuring	Range (mg/L)	Method/Chemistry
Chlorine, Free and Total	Used as disinfectant. Indicates whether residual chlorine is available to maintain proper disinfection.	0-2.00 Free 0-2.00 Total	Colorimetric/DPD
Coliform, Total and <i>E.coli</i>	Indicator organisms that can indicate possible problems with disinfection	Presence or Absence	P/A Broth with MUG
Nitrate	Can indicate non-point source pollution, the breakdown of vegetation, and oxidation of nitrogen compounds in effluents. Can be hazardous to people at levels greater than 10 mg/L (ppm).	0-30.0 mg/L	Colorimetric/ Cadmium Reduction
Nitrite	Intermediate product in the nitrogen cycle and harmful to aquatic life. Can originate from chemical fertilizers and can indicate non-point source pollution.	0-0.350 mg/L	Colorimetric/ Diazotization
Nitrogen, Ammonia	Occurs when nitrogenous products decompose in water and is harmful to aquatic life. Can indicate sewage contamination.	0-1.00 mg/L	Colorimetric/ Salicylate
pH	Common water quality indicator. Proper pH levels ensure that chlorine will perform optimally.	6.5-8.5 pH units (DR/850)	Colorimetric/ Phenol Red
		0-14 pH units	pH Pocket Pal Tester
Phosphorus, Reactive	Phosphates are used in food processing, detergents, and fertilizers. Although not directly toxic, phosphates can indicate eutrophication of water systems.	0-2.50 mg/L	Colorimetric/ Ascorbic Acid
Sulfide	Toxic to aquatic life. Can indicate industrial effluent and/or sewage contamination.	0-0.70 mg/L	Methylene Blue
TDS/ Conductivity	Can indicate high levels of salinity, metals, and inorganic compounds.	10-1999 TDS	TDS Pocket Pal Tester

1.4 Preparation for Use

Remove the instruments and accessories from the shipping boxes and inspect them for damage that may have occurred due to rough handling or extreme weather conditions. Use the Packaging Guide to verify that all the components are present. The legend is a materials list for the contents of the lab. See *Section 1.5* if any items are damaged or missing.

1.5 Missing or Damaged Product

If any items are missing or damaged, please contact the Customer Service Department in Loveland Colorado. International customers should contact the Hach office or authorized dealer serving your area. See page 93 for information on contacting Hach Company.

Please do not return items without prior authorization from Customer Service!

SECTION 2 PORTABLE INCUBATOR

Specifications are subject to change without notice.

2.1 Specifications

Ambient Operating Temperature: 0 to 40 °C

Storage Temperature: -40 to 60 °C (instrument only)

Temperature Stability: ± 0.5 °C

Temperature Range: 5° above ambient to 50 °C

Warm-up Time: 2 ± 1 hour

Capacity:

42 50-mm petri dishes

or 39 MPN tubes (19 mm OD)

or 6 P/A Disposable Bottles

Power Requirements: 12 Vdc battery or optional battery eliminator (either 115 Vac or 230 Vac available)

External Dimensions: 26 x 24 x 21 cm (10.2 x 9.4 x 8.3 in.)

Internal Dimensions: 19 x 12.5 x 13 cm (7.5 x 5 x 5.1 in.)

Instrument Weight: 1.8 kg (4 lb)

2.2 Power Selection

2.2.1 Using In A Motorized Vehicle

Plug the power cord of the instrument into the cigarette lighter outlet. The incubator will then be powered by the 12 Vdc vehicle battery.

***Note:** When using the incubator for extended periods, the automobile engine should be run periodically to ensure the automobile battery is recharged.*

2.2.2 Using Portable Battery Pack

Plug the power cord of the instrument into the battery pack. The battery pack will operate the instrument for at least 12 hours, depending on ambient temperature. To recharge the battery, plug the male plug of the recharger into a 115 Vac outlet or use the appropriate adapter for other voltage sources. The battery will completely recharge in 24 hours.

SECTION 2, continued

2.2.3 Using Battery Eliminator

Plug the power cord of the instrument into the Battery Eliminator.
Plug the battery eliminator into the AC outlet.

CAUTION

Before connecting the instrument power cord to any power source, ensure that the appropriate power supply is either 12 volt direct current (Vdc) or converts power to 12 Vdc. Before connecting the battery eliminator to a power source, ensure the appropriate line voltage is being used.

ATTENTION

Avant de raccorder le câble d'alimentation de l'appareil à toute source d'alimentation, vérifier que l'alimentation est soit une alimentation à courant continu 12 volts ou convertit l'alimentation en 12 volts continus. Avant de raccorder le transformateur d'alimentation basse tension à une source d'alimentation, s'assurer que la tension du secteur est appropriée.

PRECAUCIÓN

Antes de conectar el cable del aparato a la fuente de alimentación, asegúrese de que el suministro de corriente sea de 12 voltios de corriente continua o que procede de un transformador a 12 V CC. Antes de conectar el eliminador de baterías a la fuente de alimentación, asegúrese de que el voltaje es el adecuado.

CUATELA

Antes de ligar o cabo de alimentação a qualquer fonte de alimentação, certifique-se que a fonte de alimentação apropriada é 12 Volt corrente contínua (VCC) ou converte corrente em 12 VCC. Antes de ligar o "eliminador" de bateria a uma fonte de alimentação, assegure-se que a linha de voltagem apropriada está a ser usada.

VORSICHT


Bevor das Netzstromkabel mit einer Stromversorgung verbunden wird, überprüfen Sie bitte, ob es sich bei der jeweiligen Stromquelle entweder um 12 Volt Gleichstrom (VDC) handelt oder dass der Strom in 12 VDC umwandelt wird. Bevor man das Netzgerät an eine Stromquelle anschlieszt zu überprüfen, ob die Betriebsspannung korrekt ist.

2.3 Operation of Controls and Indicators

Figure 1 shows the Portable Incubator controls and indicators.

Key	Description
I/O	Power switch to turn instrument on and off. The switch must be on before any systems are operational, including the control circuitry.
^	Up Arrow. Press and hold until the display begins to blink, displaying the set point. Press this key to increase the incubation temperature set point.
v	Down Arrow. Press and hold until the display begins to blink, displaying the set point. Press this key to decrease the incubation temperature set point.

Figure 1 Portable Incubator Controls and Indicators



2.3.1 The Display

The display shows the current temperature inside the incubator in degrees Celsius (°C). When either the up or down arrow key is pressed and held for several seconds, the display will begin to blink, displaying the temperature set point. To change the temperature set point, press the appropriate key to increase or decrease the set point. If one of the arrow keys is not pressed within 5 seconds, the display will stop blinking and will display the current temperature inside the unit.

To re-calibrate the digital display with your reference thermometer, follow the calibration instruction in Section 2.5. The incubator was calibrated at the factory to heat to 37 °C.

2.4 Operation

CAUTION

If power sources other than those provided by Hach are used, they must be energy-limited to 12 Vdc, 15 amps maximum.

PRECAUCIÓN

Si utiliza fuentes de alimentación distintas a las suministrados por Hach, la potencia de salida deberá estar limitada a un máximo de 12 voltios, 15 amperios.

CUATELA

Se estiver a utilizar fontes de alimentação diferentes fornecidas pela Hach, estas têm de ser limitadas a um máximo de corrente de saída de 12 VCC, 15 amp.

ATTENTION

Si des sources d'alimentation autres que celles fournies par Hach sont utilisées, elles doivent être limitées à une puissance maximale de 12 volts continus, 15 ampères.

1. Plug the incubator power cord into an appropriate power source.
2. Turn the power switch on.
3. Press either the Up or Down arrow key until the display begins to blink. Use the Up or Down arrow key to adjust the display to the desired temperature set point. If one of the arrow keys is not pressed within 5 seconds, the display will stop blinking and will display the current temperature inside the unit.
4. Allow the incubator to warm up to the set point temperature.

Note: *Unit is stable when it remains at a constant temperature for 60 minutes.*

5. Open the incubator lid. Load the unit with the desired sample rack and close the lid.

2.5 Calibration

1. Plug the incubator power cord into an appropriate power source.
2. Place a reference thermometer in the center of the incubator.
3. Turn the power switch on.
4. Press either the Up or Down arrow key until the display begins to blink. Use the Up or Down arrow key to adjust the display to the desired temperature set point. If one of the arrow keys is not pressed within 5 seconds, the display will stop blinking and will display the current temperature inside the unit.
5. Allow the incubator to warm up to the set point temperature.

Note: *Unit is stable when it remains at a constant temperature for 60 minutes.*

6. Compare the reading on the reference thermometer with the digital display. If there is a difference, put the display into the calibrate mode by pressing and holding both arrow keys for about 5 seconds or until the two outside decimal points start flashing.
7. When the decimal points are flashing, use the Up or Down arrow key to adjust the display reading to match the thermometer reading.
8. Allow the unit to stabilize again.
9. Repeat the process, if necessary.

2.6 Maintenance

DANGER

For protection against electrical shock, disconnect all external electrical connections.

PERIGO

Para assegurar a protecção contra chagues eléctricos, desligue todas as ligações eléctricas externas.

PELIGRO

Como medida de protección frente a una posible descarga eléctrica, desconecte las conexiones eléctricas externa.

DANGER

Pour protection contre les chocs électriques, débrancher toutes les connexions électriques externes.

GEFAHR

Um Stromschläge zu vermeiden, müssen all elektrischen Verbindungen zum Gerät unterbrochen werden.

2.6.1 Cleaning

- Keep the incubator and accessories as clean as possible and store the instrument in the carrying case. Wipe the outside of the incubator with a soft damp cloth.
- Wipe spills up promptly.
- Clean the inside chamber with a mild soap and water solution. Rinse with clean water and wipe dry with a soft cloth.
- Foreign materials inside the unit may rust or leave spots. If corrosion is seen, scrub out the stains with a mild abrasive. Do NOT use steel wool. Failure to remove corrosion may permanently damage the liner.

2.6.2 Replacing the Fuses

The Portable Incubator itself has no fuses, but the alternate power sources available from Hach have fuses. The Battery Pack fuse is an in-line fuse and the Battery Eliminator fuses are located under the On/Off switch on the back of the Battery Eliminator enclosure.

SECTION 2, continued

2.6.2.1 Replacing the 115/230 Volt Battery Eliminator Fuses

This applies to Hach Cat. No. 25804-00

1. Disconnect all power cords from the Battery Eliminator.
2. Use a small flat screwdriver or blade to gently pry the faceplate out of the enclosure. Pull on the faceplate until the fuse holder comes out.
3. The fuse holder has two T, 2 Amp, 250 V fuses. Both fuses must be operational for the Battery Eliminator to work. Remove blown fuses by gently pulling them out of the holder. Replace with a good fuse.
4. Replace the fuse holder into the instrument.

2.6.2.2 Replacing the Battery Pack Fuse

1. Disconnect the Battery Pack cord from the Battery Pack.
2. The in-line fuse is located 2 to 3 inches from the Battery Pack connection in a plastic holder.
3. Twist the upper and lower part of the fuse holder in opposite directions to open the fuse holder.
4. Remove the ceramic fuse and replace it with a 15 Amp, 250 V fuse.
5. Reverse Steps 1 through 3 to reassemble.

REPLACEMENT PARTS

Description	Unit	Cat. No.
Fuse, T, 2A, 250 V	1	45686-00
Fuse, 15 A, 250 V, ceramic	1	25850-00
Hinge	1	25848-00
Incubator, portable, for MELs	1	25699-00
Lid Assembly	1	25849-00

PORTABLE INCUBATOR ACCESSORIES

Battery Eliminator, 115/230 Vac (unlisted).....	1	25804-00
Battery Eliminator, 230 Vac, CE/VDE approved	1	40277-00
Battery, portable, 12 Vdc, rechargeable	1	25803-01
Battery, portable, 12 Vdc rechargeable, with carrying case and 115 Vac recharger.....	1	25803-00
Battery Recharger Adapter, 12 Vdc, 230 Vac, CE/VDE, for 25803-01	1	25959-01
Battery Recharger Adapter, 12 Vdc, 230 Vac, European plug, unlisted	1	25959-02

SECTION 2, continued

PORTABLE INCUBATOR ACCESSORIES, continued

Description	Unit	Cat. No.
Rack, MF/General Purpose, holds 42 50-mm petri dishes	1	25805-02
Rack, MPN tube, holds 39 tubes, 19-mm OD	1	25805-01
Rack, P/A bottle, holds 6 bottles, 5-cm OD.....	1.....	25805-00
Sample Transport Kit: includes cooler, plastic-coated rack, 100 sampling bags with dechlorinating agent, and a refrigerant pack	1.....	25687-00

MEDIA SETS (includes media and consumables required for testing)

MEL/MF Media Set for Total Coliforms (200 tests)

Includes:

200 m-Endo Broth PourRite Ampules, 200 sterile 50-mm petri dishes with pads,
216 sterile 0.45- μ m membrane filters, 216 sterile push-fit funnels,
200 Whirl-Pak bags with dechlorinating agent 25801-00

MEL/MPN Media Set for Total Coliforms and *E. coli* (25 5-tube tests)

Includes:

135 LT/MUG Broth tubes, 30 BGB Broth tubes (for total coliform confirmation),
50 sterile inoculating loops, 25 Whirl-Pak bags with dechlorinating agent) 25802-00

MEL P/A Media Set for Total Coliforms and *E. coli*

See “REQUIRED MEDIA AND APPARATUS” on page 24.

MEL/850 Portable Water Laboratory

MEL/850 Reagent Set

Includes:

Reagents for measuring Ammonia Nitrogen, Coliform Bacteria, Free Chlorine,
Total Chlorine, Nitrate, Nitrite, pH, Phosphate, and Sulfide (100 tests each) 26911-00

MEL/850 Apparatus Set

Includes:

Clippers, pocket thermometer, demineralizer bottle, 2 100-mL beakers,
portable UV lamp, P/A bottle rack, TDS Pocket Pal Tester, pH Pocket Pal
Tester, 25-mL graduated plastic cylinder, 100 Whirl-Pak bags with
dechlorinating agent, 60 plastic droppers with 0.5- and 1.0-mL marks..... 26912-00

For Technical Assistance, Price and Ordering

In the U.S.A.—Call 800-227-4224

Outside the U.S.A.—Contact the Hach office or distributor serving you.

3.1 Overview of Coliform Bacteria

Many of the microorganisms that cause serious disease, such as typhoid fever, cholera, and dysentery, can be traced directly to polluted water. These disease-producing organisms, or pathogens, are discharged along with fecal wastes and are difficult to detect in water supplies. People may contact these pathogens in swimming pools, on bathing beaches, in rivers and streams, and from drinking contaminated water.

Testing for bacterial pathogens in water is impractical for a number of reasons, such as lengthy and involved test procedures. Most microbiological testing of water measures indicator organisms, not pathogens. Indicator organisms are bacteria that may not be pathogenic but usually are present when pathogens are present, and are more resistant to environmental stresses than pathogens. No organism or group of organisms satisfies all of the criteria for an indicator; however, coliforms satisfy most of the requirements.

Total coliform tests are used for potable water supplies. Fecal coliform tests usually are performed on untreated non-potable water, wastewater, bathing water, and swimming water.

For simultaneous detection of total coliforms and *Escherichia coli* (*E. coli*), a type of fecal coliform, Hach offers m-ColiBlue24[®] Broth, Presence/Absence Bromcresol Purple Broth with MUG, and Lauryl Tryptose with MUG Broth.

The Presence/Absence (P/A) method is a qualitative test that indicates only the presence or absence of organisms, not the number of organisms. The P/A method is fast and suited to spot-checking applications. Only a minimal amount of analytical experience is required to perform the test. Simply combine sample with medium, incubate for 24 to 48 hours, and check for a reaction indicating the presence of either total coliforms or *E. coli*.

U.S. Environmental Protection Agency drinking water regulations, effective December 31, 1990, require reporting only the presence or absence of coliforms. The World Health Organization recommends using the P/A method for drinking water to ensure zero total coliforms and zero fecal coliforms or *E. coli*. The maximum contaminant goal of zero total coliforms eliminates the need to enumerate coliforms.

3.2 Testing Techniques

Good laboratory technique is essential when accuracy is important, particularly in microbiological laboratory procedures. To assure reliable results use care when collecting and preserving samples; carefully clean the laboratory or work surface; use proper sterilization and inoculation practices; and maintain close temperature control.

Using high-quality laboratory equipment and ready-to-use media can also save time and minimize errors. Hach's prepared media help eliminate contamination due to individual technician technique.

3.3 Preparing Sample Containers

Take care to prevent contamination when conducting bacterial tests. All materials used for containing or transferring samples must be sterile! To collect samples, use presterilized plastic bags, presterilized disposable bottles, autoclavable glass bottles, or autoclavable plastic bottles.

3.3.1 Presterilized Containers

Presterilized plastic bags with dechlorinating agent are included in the MEL and media sets. Plastic bags are available presterilized with or without dechlorinating agent (sodium thiosulfate). Presterilized bottles are available with a 100-mL fill-to line.

Note: Use dechlorinating agent with potable water or chlorinated water samples. It is not necessary for unchlorinated water samples. However, dechlorinating agent will not interfere with unchlorinated samples. For simplicity, plastic bags containing dechlorinating agent may be used for all samples.

3.3.2 Autoclavable Containers

Glass or plastic bottles (at least 125 mL) may be used instead of plastic bags or bottles. Prepare these containers as follows:

1. Wash in hot water and detergent.
2. Thoroughly rinse with hot tap water, followed by a distilled water rinse to make sure that all detergent is removed.
3. If dechlorinating agent is needed (for chlorinated, potable water) add one Dechlorinating Reagent Powder Pillow to each 125-mL sample container. Add two powder pillows to a 250-mL sample container.
4. Steam sterilize glass and autoclavable plastic containers at 121 °C for 15 minutes. Glass sample containers may be sterilized by hot air at 170 °C for 1 hour. Store sterile containers tightly capped in a clean environment until needed.

3.4 Collecting and Preserving Samples

Proper sampling technique ensures that seasonal variances are detected and that results are representative of the sample source.

Collect a sufficient volume of water for analysis, at least 100 mL of sample. World Health Organization guidelines suggest 200 mL per sample, while *Standard Methods for the Examination of Water and Wastewater* guidelines suggest 100 mL per sample. Avoid sample contamination during collection.

Dechlorination is not necessary if the sample is added directly to the medium on site. Otherwise, samples should be dechlorinated and immediately transported for analysis. Sodium thiosulfate, which has been sterilized within the collection vessel, is generally used to destroy chlorine residual.

Analyze as soon as possible after collection. The maximum time between collection and examination of samples should be 8 hours (maximum transit time 6 hours, maximum processing time 2 hours). *If the time between collection and analysis will exceed 8 hours, maintain the sample at/or below 4 °C, but do not freeze.* Maximum time between collection and analysis should not exceed 24 hours. Failure to properly collect and transport samples will cause inaccurate results.

Potable water should contain no coliforms per 100 mL, so testing should be done on undiluted samples. Collect at least 100 mL of sample in presterilized plastic bags or bottles or in sterile glass or plastic sample bottles. Sample containers should not be filled completely. Maintain at least 2.5 cm (approximately 1 in.) of air space to allow adequate space for mixing the sample prior to analysis.

MEL Portable Incubator Laboratories allow you to conduct all of your analysis in the field. But, if you are just transporting samples to a lab, your samples should arrive at the lab within 24 hours after collection. In warm climates, the samples must be packed in a freezing mixture to maintain the sample temperature between 4 and 10 °C. Failure to properly collect and transport samples will cause inaccurate results.

SECTION 3, continued

3.4.1 Samples from Faucets, Spigots, Hydrants, or Pumps

Collect representative samples by allowing water to run from a faucet, spigot, hydrant, or pump at a moderate rate (without splashing) for 2 to 3 minutes before sampling. Do not adjust the flow rate while collecting the sample. Valves, spigots and faucets that swivel or leak should be avoided. Remove attachments, such as aerators and screens, prior to sampling.

Carefully open sample containers just prior to collection and close immediately following collection. Do not lay the lid or cap down. Avoid touching the mouths and insides of the containers. Do not rinse the containers. Properly label each sample container and analyze samples as soon as possible after collection.

3.4.2 Samples from Rivers, Lakes, and Reservoirs

When sampling a river, lake, or reservoir, fill the sample container below the water surface. Do not sample near the edge or bank. Remove the cap, grasp the sample container near the bottom and plunge the container, mouth down, into the water. (This technique excludes any surface scum.) Fill the container by positioning the mouth into the current or, in nonflowing water, by tilting the bottle slightly and allowing it to fill slowly. Do not rinse. Label each sample container and analyze samples as soon as possible after collection.

3.5 Disposing of Completed Tests

Active bacterial cultures grown during incubation must be disposed of safely. This may be accomplished in one of two ways.

1. **Bleach.** Used test containers may be sterilized by using a 10% bleach solution. Pour the test container contents and the test containers into the bleach solution. Allow 10 to 15 minutes contact time with the bleach. Pour the liquid down the drain. Dispose of the test containers in the normal waste.
2. **Autoclave.** Place used test containers in a contaminated items bag or a biohazard bag and seal tightly. Test containers must be placed in a bag before autoclaving to prevent media leakage into the autoclave. Autoclave test containers at 121 °C for 15 minutes at 15 pounds pressure. Once sterile, dispose of the test containers in the normal waste. Place the bag of test containers in a separate garbage bag and tie tightly.

3.6 Presence/Absence Methods—Bromcresol Purple Broth (Method 8319) and Bromcresol Purple Broth with MUG (Method 8364)

3.6.1 Preparing Materials

To save time, start the incubator while preparing other materials. Set the incubator for the proper temperature setting described in the procedure (usually total coliforms are incubated at 35 ± 0.5 °C and fecal coliforms are incubated at 45 ± 0.2 °C).

Disinfect the work bench with a germicidal cloth, dilute bleach solutions, bactericidal spray, or dilute iodine solution. Wash hands thoroughly with soap and water.

Mark container with the sample number, dilution, date, and other necessary information. Take care not to contaminate the inside of the sample container in any way.

See the preceding pages for information about preparing sample containers, and collecting and preserving samples.

3.6.2 Convenient Packaging

Hach's Presence/Absence (P/A) Broth and P/A Broth with MUG come packaged in disposable bottles and in glass ampules. To reduce medium costs, both bottles and ampules contain 20 mL of 6X strength medium. The medium is sterilized by membrane filtration to prevent degradation. Disposable bottles and ampules are shipped with a Certificate of Analysis and have an expiration or manufacture date printed on the label.


SECTION 3, continued

3.7 Using P/A Broth Disposable Bottles




1. Collect 100 mL of sample in a sterile container. Do not contaminate the sample or sample container.

Note: Remove screens and other aeration devices from faucets and let water run for 2 to 3 minutes before collecting the sample.



2. Add sample to the fill line on the P/A Broth Sample Bottle. Sample may be added from a sterile container, or directly from a faucet or spigot.




3. Incubate the samples at 35 ± 0.5 °C for 24 to 48 hours.



4. Note the reaction after 24 hours of incubation. If sample is negative, continue incubating for another 24 hours. See *Table 1*.



5. Confirm positive samples by inoculating the appropriate media from positive P/A Broth samples. See *Table 2*.




6. Dispose of completed tests appropriately. See *Disposing of Completed Tests* on page 18.

3.8 Using P/A Broth Ampules




1. Collect 100 mL of sample in a sterile container. Do not contaminate the sample or sample container.


Note: Remove screens and other aeration devices from faucets and let water run for 2 to 3 minutes before collecting the sample.



2. Add sample to the fill line on the P/A Broth Sample Bottle. Sample may be added from a sterile container, or directly from a faucet or spigot.



3. Add the contents of one P/A Broth Ampule to the 100 mL of sample.




4. Incubate the samples at 35 ± 0.5 °C for 24 to 48 hours.



5. Note the reaction after 24 hours of incubation. If sample is negative, continue incubating for another 24 hours. See *Table 1*.



6. Confirm positive samples by inoculating the appropriate media from positive P/A Broth samples. See *Table 2*.



7. Dispose of completed tests appropriately. See *Disposing of Completed Tests* on page 18.

SECTION 3, continued

3.9 Interpreting P/A Results

Table 1 Reactions Using P/A Broth

Reaction	Comments	Report as:
Color change from reddish purple to yellow or yellow brown		Positive for total coliforms
No color change after 24 hours	Incubate for another 24 hours and re-check the sample for color change	
No color change after 48 ± 3 hours		Negative for total coliforms
Fluorescence under long-wave UV light (if using P/A Broth with MUG)		Positive for <i>E. coli</i>

3.10 Confirming Positive Samples

Inoculum from incubated samples can be used to confirm the presence of bacteria. See *Table 2*. The media listed for fecal coliforms, total coliforms, and *E. coli* are USEPA-accepted for reporting purposes.

Table 2 Confirmation Media

Bacteria	Confirmation Media	Incubation	Positive Result
Total Coliforms (USEPA)	Brilliant Green Bile Broth (Cat. No. 322-15; 15/pkg)	24-48 hours at 35 ± 0.5 °C	Gas
Fecal Coliforms (USEPA)	EC Medium (Cat. No. 14104-15; 15/pkg)	24 hours at 44.5 ± 0.2 °C	Gas
<i>E. coli</i> (USEPA)	EC Medium with MUG (Cat. No. 24715-15; 15/pkg)	24 hours at 44.5 ± 0.2 °C	Gas = positive for fecal coliforms Fluorescence = positive for <i>E. coli</i>

3.11 Summary of Method 8319

Hach's Presence/Absence (P/A) Bromcresol Purple Broth is ideal for screening drinking water samples for total coliforms. The method is a simple modification of the multiple-tube method. It uses lactose and lauryl tryptose broths with bromcresol purple, which detects acidity formed during lactose fermentation by the bacteria.

Simply combine 100 mL of sample and P/A Broth, incubate for 24 hours and check for a color change. A yellow or yellow-brown color indicates the presence of total coliforms.

3.12 Summary of Method 8364

Hach's P/A Bromcresol Purple Broth with MUG allows simultaneous detection of total coliform bacteria and *E. coli*. In addition to the lactose and lauryl tryptose broths with bromcresol purple, this medium contains MUG reagent (4-methylumbelliferyl- β -D-glucuronide). MUG reagent produces a fluorogenic product when hydrolyzed by glucuronidase (an enzyme specific to *E. coli*). MUG detects non-gas producing (anaerogenic) strains of *E. coli* and works well when competitive organisms are present.

Simply combine 100 mL of sample and P/A Broth with MUG, incubate for 24 hours and check for a color change and fluorescence. A yellow or yellow-brown color indicates the presence of total coliforms. To detect *E. coli*, examine samples under a long-wave ultraviolet (UV) light. Fluorescence indicates the presence of *E. coli*.

SECTION 3, continued

REQUIRED MEDIA AND APPARATUS

Description	Quantity Required		Unit	Cat No.
	Per Test			
Bags, Sampling, Sterile Whirl-Pak with dechlorinating agent, 170-mL.....	1	100/pkg.....	20753-33
Incubator, portable, for MELs	1	each.....	25699-00
Presence/Absence Broth w/MUG, disposable bottles ..	1	50/pkg.....	24016-50
P/A Bottle Rack (for use with Portable Incubator).....	1	each.....	25805-00
UV Lamp, long-wave, portable, 4 watt.....	1	each.....	24152-00

OPTIONAL MEDIA AND APPARATUS

Bags, for contaminated items.....	200/pkg.....	24633-00
Bottles, presterilized, 100-mL fill-to line	12/pkg.....	24950-12
Bottles, presterilized, 100-mL fill-to line	50/pkg.....	24950-50
Breaker, P/A Ampule.....	1 25640-00
Germicidal Cloths.....	50/pkg.....	24632-00
Inoculating Loops, sterile 10 µL disposable (for confirmation)	25/pkg.....	27491-25
P/A Broth Ampules.....	25/pkg.....	24949-25
P/A Broth Ampules with MUG	25/pkg.....	24955-25
P/A Broth Disposable Bottles	12/pkg.....	23232-12
P/A Broth Disposable Bottles	50/pkg.....	23232-50
P/A Broth Disposable Bottles with MUG.....	12/pkg.....	24016-12
Pen, Laboratory.....	1 20920-00
UV Lamp, long-wave, 115 Vac.....	1 21843-00
UV Lamp, long-wave, 230 Vac.....	1 21843-02

For Technical Assistance, Price and Ordering

In the U.S.A.—Call 800-227-4224

Outside the U.S.A.—Contact the Hach office or distributor serving you.

SECTION 4 MEASURING TOTAL DISSOLVED SOLIDS (TDS)

4.1 Calibrating the TDS Pocket Pal Tester

Verify the accuracy of the tester before use and periodically thereafter as follows:

1. Add a solution with a known TDS value to a 50-mL beaker.

Note: Best accuracy ($\pm 2\%$) is obtained when samples are measured at the same temperature ($25\text{ }^{\circ}\text{C}$) as that of the standards used for calibration. If samples are measured at a different temperature, the Pocket Pal Tester compensates for the difference by adjusting the reading $2\% \text{ }^{\circ}\text{C}$. The accuracy of temperature-compensated readings is $\pm 10\%$ in the temperature range of 0 to $50\text{ }^{\circ}\text{C}$.

2. Press the ON/OFF switch once to turn the tester on. (See Figure 2.)
3. Remove the protective cap from the bottom.
4. Immerse the bottom of the tester 2.5 to 8.9 cm (1 to 3.5 in.) into the standard.
5. Using the tester, gently stir the standard for several seconds. When the digital display stabilizes, read the TDS value.

Note: Readings may not stabilize for up to 2 minutes, especially if the temperature is far from ambient.

6. If necessary, adjust the Calibration Trimmer (see Figure 3) using the supplied trimmer tool (or a small flat-bladed screwdriver). Turn the trimmer tool until the reading corresponds to TDS value of the standard.

Note: To maintain or improve performance, periodically clean the stainless steel electrodes by rinsing in isopropyl alcohol.

Figure 2






Figure 3



4.2 Replacing the Batteries

1. Use a coin to turn the battery compartment cover, located on the top of the tester, to the left 1/4 turn. (see *Figure 4*).
2. Remove the cover.
3. Replace all four batteries with EverReady E675E, Duracell RM675, or Hach batteries, Cat. No. 23678-00, in the same orientation (polarity) as they were removed.
4. Replace the cover.

Figure 4




SECTION 4, continued


4.3 TDS Procedure



1. Pour sample into a 50-mL plastic beaker.




2. Immerse the bottom of the tester 2.5 to 8.9 cm (1.0 to 3.5 in.) into the sample. Stir gently for several seconds.




3. When the digital display stabilizes, read and record the TDS value.

Note: Readings may not stabilize for up to 2 minutes, especially if the temperature is far from ambient.



4. Rinse the bottom of the tester with deionized water. Wipe with a tissue before continuing to the next sample.



5. After testing, rinse and dry the tester, press the on/off switch to off, and replace the cap.

Note: To maintain or improve performance, periodically clean the stainless steel electrodes by rinsing in isopropyl alcohol.

SECTION 4, continued

REQUIRED REAGENTS AND APPARATUS

Description	Quantity Required		Unit	Cat. No.
	Per Test			
Batteries, 1.4 V (offer 100 hours of continuous use)....	4	4/pkg 23678-00
Beaker, polypropylene, 50 mL.....	1	each 1080-41
TDS Pocket Pal Tester, 10 to 1990 TDS.....	1	each 44400-01
TDS Standard, sodium chloride, 85.47 mg/L NaCl..	varies	100 mL 23075-42
TDS Standard, sodium chloride, 491 mg/L NaCl.....	varies	100 mL 14400-42
TDS Standard, sodium chloride, 1000 mg/L NaCl...	varies	100 mL 2105-42

For Technical Assistance, Price and Ordering

In the U.S.A.—Call 800-227-4224

Outside the U.S.A.—Contact the Hach office or distributor serving you.

SECTION 5 MEASURING pH

5.1 Calibrating the pH Pocket Pal Tester

Verify the accuracy of the tester before use and periodically thereafter as follows:

1. Prepare a pH 7.00 and a pH 4.01 or 10.00 buffer (Cat. No. 22270-66, 22269-66, or 22271-66, respectively). Add a Buffer Powder Pillow to 50 mL of water. Mix well.

Note: Best accuracy is obtained when samples are measured at the same temperature (25 °C) as that of the standards used for calibration.

2. Press the ON/OFF switch once to turn the tester on. See *Figure 5*.
3. Remove the protective cap from the bottom.
4. Immerse the bottom of the tester 2.5 to 8.9 cm (1 to 3.5 in.) into the standard.
5. Using the tester, gently stir the standard for several seconds. When the digital display stabilizes, read the pH value.
6. If necessary, adjust the Calibration Trimmer (see *Figure 6*) using the supplied trimmer tool (or a small flat-bladed screwdriver). Turn the trimmer tool until the reading corresponds to the pH value of the buffer (7.0 or 4.0/10.0 pH).
7. Rinse the bottom of the tester. Replace the cap.

Figure 5






Figure 6



5.2 Replacing the Batteries

1. Use a coin to turn the battery compartment cover, located on the top of the tester, to the left 1/4 turn. (see *Figure 7*).
2. Remove the cover.
3. Replace all four batteries with EverReady E675E, Duracell RM675, or Hach batteries, Cat. No. 23678-00, in the same orientation (polarity) as they were removed.
4. Replace the cover.


Figure 7




5.3 pH Procedure



1. Pour sample into a 50-mL plastic beaker.




2. Immerse the bottom of the tester 2.5 to 8.9 cm (1.0 to 3.5 in.) into the sample. Stir gently for several seconds.




3. When the digital display stabilizes, read and record the pH value.

Note: If pH readings are erratic, replace the batteries (see Section 5.2).



4. Rinse the bottom of the tester with deionized water. Wipe with a tissue before continuing to the next sample.



5. After testing, rinse and dry the tester, press the on/off switch to off, and replace the cap.

Note: For faster response and longer tester life, place deionized water in the cap to prevent the glass bulb from drying out.

Note: Soak the electrode tip in tap water for a few minutes each week to condition the electrode.

Note: Potassium chloride, used as reference solution electrolyte, may deposit on the tester as white precipitate. This does not affect performance; remove it with a damp cloth or tissue.

SECTION 5, continued

REQUIRED REAGENTS AND APPARATUS

Description	Quantity Required		Unit	Cat. No.
	Per Test			
Batteries, 1.4 V (offer 100 hours of continuous use)....	4	4/pkg.....	23678-00
Beaker, polypropylene, 50 mL.....	1	each.....	1080-41
Buffer Powder Pillows, pH 4.01	1	50/pkg.....	22269-66
Buffer Powder Pillows, pH 7.00	1	50/pkg.....	22270-66
Buffer Powder Pillows, pH 10.00	1	50/pkg.....	22271-66
pH Pocket Pal Tester, 0 to 14 pH units	1	each.....	44350-01

For Technical Assistance, Price and Ordering

In the U.S.A.—Call 800-227-4224

Outside the U.S.A.—Contact the Hach office or distributor serving you.




DR/850 COLORIMETER PROCEDURES

CHLORINE, FREE (0 to 2.00 mg/L)

For water, wastewater, and seawater

DPD Method (Powder Pillows or AccuVac Ampuls)

USEPA accepted for reporting wastewater and drinking water analyses*

Using Powder Pillows


1. Enter the stored program number for free and total chlorine (Cl₂) powder pillows.

Press: **PRGM**

The display will show:

PRGM ?

Note: For most accurate results, perform a Reagent Blank Correction using deionized water (see Section 1 of the DR/850 Procedures Manual).




2. Press: **9 ENTER**

The display will show **mg/L, Cl₂** and the **ZERO** icon.



3. Fill a sample cell with 10 mL of sample (the blank).


Note: Samples must be analyzed immediately and cannot be preserved for later analysis.



4. Place the blank into the cell holder. Tightly cover the sample cell with the instrument cap.

* Procedure is equivalent to USEPA method 330.5 for wastewater and Standard Method 4500-Cl G for drinking water.

CHLORINE, FREE, continued



5. Press: **ZERO**


The cursor will move to the right, then the display will show:

0.00 mg/L Cl₂

Note: If Reagent Blank Correction is on, the display may flash "limit". See Section 1 of the DR/850 Procedures Manual.




6. Fill another cell with 10 mL of sample.




7. Add the contents of one DPD Free Chlorine Powder Pillow to the sample cell (the prepared sample). Cap the cell and swirl vigorously to dissolve the powder.

Note: A pink color will develop if free chlorine is present.



8. Immediately place the prepared sample into the cell holder. Tightly cover the sample cell with the instrument cap.

Note: Perform Step 9 within 1 minute of reagent addition.




9. Press: **READ**

The cursor will move to the right, then the result in mg/L chlorine will be displayed.

Note: Standard Adjust may be performed using a prepared standard (see Section 1 of the DR/850 Procedures Manual).

Note: If the sample temporarily turns yellow after reagent addition, or the display flashes "limit", it is due to high chlorine levels. Dilute a fresh sample and repeat the test. A slight loss of chlorine may occur during dilution. Multiply the result by the dilution factor; see Section 1 of the DR/850 Procedures Manual.

Using AccuVac Ampuls




1. Enter the stored program number for free and total chlorine (Cl_2) AccuVac Ampuls.

Press: **PRGM**

The display will show:


PRGM ?

Note: For most accurate results, perform a Reagent Blank Correction using deionized water (see Section 1 of the DR/850 Procedures Manual).




2. Press: **11 ENTER**

The display will show **mg/L, Cl2** and the **ZERO** icon.




3. Fill a sample cell with at least 10 mL of sample (the blank). Collect at least 40 mL of sample in a 50-mL beaker.

Note: Samples must be analyzed immediately and cannot be preserved for later analysis.



4. Place the blank into the cell holder. Tightly cover the sample cell with the instrument cap.

CHLORINE, FREE, continued




5. Press: **ZERO**

The cursor will move to the right, then the display will show:


0.00 mg/L Cl₂

Note: If Reagent Blank Correction is on, the display may flash "limit". See Section 1 of the DR/850 Procedures Manual.




6. Fill a DPD Free Chlorine Reagent AccuVac Ampul with sample.

Note: Keep the tip immersed while the ampul fills completely.




7. Quickly invert the ampul several times to mix. Wipe off any liquid or fingerprints.

Note: A pink color will form if chlorine is present.



8. Immediately place the AccuVac Ampul into the cell holder. Tightly cover the ampul with the instrument cap.

Note: Perform Step 9 within 1 minute of reagent addition.



9. Press: **READ**

The cursor will move to the right, then the result in mg/L chlorine will be displayed.

Note: Standard Adjust may be performed using a prepared standard (see Section 1 of the DR/850 Procedures Manual).

Note: If the sample temporarily turns yellow after reagent addition, or the display flashes "limit", it is due to high chlorine levels. Dilute a fresh sample and repeat the test. A slight loss of chlorine may occur during dilution. Multiply the result by the dilution factor; see Section 1 of the DR/850 Procedures Manual.

Sampling and Storage

Analyze samples for chlorine **immediately** after collection. Free chlorine is a strong oxidizing agent, and it is unstable in natural waters. It reacts rapidly with various inorganic compounds and more slowly oxidizes organic compounds. Many factors, including reactant concentrations, sunlight, pH, temperature, and salinity influence decomposition of free chlorine in water.

Avoid plastic containers since these may have a large chlorine demand. **Pretreat glass** sample containers to remove any chlorine demand by soaking in a dilute bleach solution (1 mL commercial bleach to 1 liter of deionized water) for at least 1 hour. Rinse thoroughly with deionized or distilled water. If sample containers are rinsed thoroughly with deionized or distilled water after use, only occasional pretreatment is necessary.

Do not use the same sample cells for free and total chlorine. If trace iodide from the total chlorine reagent is carried over into the free chlorine determination, monochloramine will interfere. It is best to use separate, dedicated sample cells for free and total chlorine determinations.

A common error in testing for chlorine is introduced when a representative sample is not obtained. If sampling from a tap, let the water flow for at least 5 minutes to ensure a representative sample. Let the container overflow with the sample several times, then cap the sample container so there is no headspace (air) above the sample. If sampling with a sample cell, rinse the cell several times with the sample, then carefully fill to the 10-mL mark. Perform the analysis immediately.

Accuracy Check

Standard Additions Method (using powder pillows)

- a) Snap the top off a LR Chlorine PourRite Ampule Standard Solution.
- b) Use a TenSette Pipet to add 0.1 mL of the standard to the reacted sample (this is the spiked sample). Swirl to mix.
- c) Re-zero the instrument using the original sample (the blank).
- d) Place the spiked sample in the cell holder and press **READ**. Record the results.

CHLORINE, FREE, continued

- e) Calculate the concentration of mg/L chlorine added to the sample:

$$\text{mg/L chlorine added} = \frac{0.1(\text{vol. standard added}) \times \text{Certificate value (mg/L Cl}_2\text{)}}{10.1(\text{sample} + \text{standard volume})}$$

- f) The spiked sample result (*Step d*) should reflect the analyzed sample result + the calculated mg/L Cl₂ added (*Step e*).
- g) If this increase does not occur, see *Standard Additions* in *Section 1* of the *DR/850 Procedures Manual* for more information.

Standard Additions Method (using AccuVac Ampuls)

- a) Snap the top off a LR Chlorine PourRite Ampule Standard Solution.
- b) Use a graduated cylinder to measure 25 mL of sample into each of two beakers.
- c) Use a TenSette Pipet to add 0.2 mL of the standard to one of the beakers (this is the spiked sample). Swirl to mix.
- d) Fill a DPD Free Chlorine AccuVac completely from each beaker.
- e) Analyze the spiked and unspiked sample as described in the procedure.
- f) Calculate the concentration of mg/L chlorine added to the sample:

$$\text{mg/L chlorine added} = \frac{0.2(\text{vol. standard added}) \times \text{Certificate value (mg/L Cl}_2\text{)}}{25.2(\text{sample} + \text{standard volume})}$$

The spiked sample result should reflect the analyzed sample result + the calculated mg/L Cl₂ added (*Step f*).

- g) If this increase does not occur, see *Standard Additions* in *Section 1* of the *DR/850 Procedure Manual* for more information.

Method Performance

Precision

In a single laboratory using a standard solution of 1.00 mg/L chlorine and two representative lots of reagents with the instrument, a single operator obtained a standard deviation of \pm 0.01 mg/L chlorine.

In a single laboratory using a standard solution of 1.00 mg/L chlorine and two representative lots of AccuVac ampuls with the instrument, a single operator obtained a standard deviation of \pm 0.01 mg/L chlorine.

Estimated Detection Limit (EDL)

The estimated detection limit for programs 9 and 11 is 0.02 mg/L Cl₂. For more information on derivation and use of Hach's estimated detection limit, see *Section 1* of the *DR/850 Procedure Manual*.

CHLORINE, FREE, continued

Interferences

Interfering Substance	Interference Level and Treatment
Acidity	Greater than 150 mg/L CaCO ₃ . May not develop full color or color may fade instantly. Neutralize to pH 6-7 with 1 N sodium hydroxide. Determine amount to be added on separate sample aliquot, then add the same amount to the sample being tested. Correct for volume addition (see <i>Correcting for Volume Additions, Section 1</i> of the <i>DR/850 Procedure Manual</i>).
Alkalinity	Greater than 250 mg/L CaCO ₃ . May not develop full color or color may fade instantly. Neutralize to pH 6-7 with 1 N sulfuric acid. Determine amount to be added on separate sample aliquot, then add the same amount to the sample being tested. Correct for volume addition (<i>Correcting for Volume Additions, Section 1</i> of the <i>DR/850 Procedure Manual</i>).
Bromine	Interferes at all levels
Chlorine Dioxide	Interferes at all levels
Chloramines, organic	May interfere
Hardness	No effect at less than 1,000 mg/L as CaCO ₃
Iodine	Interferes at all levels
Manganese, Oxidized (Mn ⁴⁺ , Mn ⁷⁺) or Chromium, Oxidized (Cr ⁶⁺)	<ol style="list-style-type: none">1. Adjust sample pH to 6-7.2. Add 3 drops potassium iodide (30 g/L) to a 25-mL sample.3. Mix and wait 1 minute.4. Add 3 drops sodium arsenite (5 g/L) and mix.5. Analyze 10 mL of the treated sample as described in the procedure.6. Subtract the result from this test from the original analysis to obtain the correct chlorine concentration.
Monochloramine	Causes a gradual drift to higher readings. When read within 1 minute after reagent addition, 3 mg/L monochloramine causes less than a 0.1 mg/L increase in the reading.
Ozone	Interferes at all levels
Peroxides	May interfere
Extreme sample pH and highly buffered samples	Adjust to pH 6-7. See <i>Interferences, Section 1</i> of the <i>DR/850 Procedure Manual</i> .

CHLORINE, FREE, continued

Summary of Method

Chlorine in the sample as hypochlorous acid or hypochlorite ion (free chlorine or free available chlorine) immediately reacts with DPD (N,N-diethyl-p-phenylenediamine) indicator to form a magenta color which is proportional to the chlorine concentration.

REQUIRED REAGENTS & APPARATUS (Using Powder Pillows)

Description	Quantity Required		Unit	Cat. No.
	Per Test			
DPD Free Chlorine Powder Pillows, 10 mL	1 pillow.....		100/pkg	21055-69
Sample Cell, 10-20-25 mL, with cap	2		6/pkg	24019-06

REQUIRED REAGENTS & APPARATUS (Using AccuVac Ampuls)

DPD Free Chlorine Reagent AccuVac Ampuls	1 ampul.....		25/pkg	25020-25
Beaker, 50 mL	1		each	500-41H

OPTIONAL REAGENTS

Chlorine Standard Solution, PourRite Ampule, 25-30 mg/L, 2mL	20/pkg			26300-20
DPD Free Chlorine Reagent, Swiftest Dispenser.....	250 tests			28023-00
Potassium Iodide Solution, 30 g/L	100 mL*	MDB		343-32
Sodium Arsenite, 5 g/L	100 mL*	MDB		1047-32
Sodium Hydroxide Standard Solution, 1.00 N	100 mL*	MDB		1045-32
Sulfuric Acid Standard Solution, 1.000 N.....	100 mL*	MDB		1270-32
Water, deionized	4L			272-56

OPTIONAL APPARATUS

AccuVac Snapper	each			24052-00
Cylinder, graduated, 25 mL	each			508-40
pH Meter, portable	each			HQ11d
pH Electrode, gel, standard, with 1 meter cable	each			PHC10101
pH Paper, 1 to 11 pH units	5 rolls/pkg			391-33
Pipet, TenSette, 0.1 to 1.0 mL	each			19700-01
Pipet Tips, for 19700-01 TenSette Pipet	50/pkg			21856-96
PourRite Ampule Breaker	each			24846-00

For Technical Assistance, Price and Ordering

In the U.S.A.—Call 800-227-4224

Outside the U.S.A.—Contact the Hach office or distributor serving you.


* Marked Dropper Bottle - contact Hach for larger sizes.

CHLORINE, TOTAL (0 to 2.00 mg/L)

For water, wastewater, and seawater

DPD Method (Powder Pillows or AccuVac Ampuls)

USEPA accepted for reporting water and wastewater analyses*

Using Powder Pillows


1. Enter the stored program number for total chlorine (Cl₂) powder pillows.

Press: **PRGM**

The display will show:


PRGM ?

Note: For most accurate results, perform a Reagent Blank Correction using deionized water (see Section 1 of the DR/850 Procedures Manual).




2. Press: **9 ENTER**

The display will show **mg/L, Cl₂** and the **ZERO** icon.



3. Fill a sample cell with 10 mL of sample (the blank).


Note: Samples must be analyzed immediately and cannot be preserved for later analysis.



4. Place the blank into the cell holder. Tightly cover the sample cell with the instrument cap.

* Procedure is equivalent to USEPA method 330.5 for wastewater and Standard Method 4500-Cl G for drinking water.

CHLORINE, TOTAL, continued



5. Press: ZERO


The cursor will move to the right, then the display will show:

0.00 mg/L Cl₂

Note: If Reagent Blank Correction is on, the display may flash "limit". See Section 1 of the DR/850 Procedures Manual.




6. Fill a second cell to the 10-mL mark with sample.



7. Add the contents of one DPD Total Chlorine Powder Pillow to the sample cell (the prepared sample). Cap and swirl the sample cell vigorously to dissolve the powder.

Note: It is not necessary that all the powder dissolves.




8. Press:

TIMER ENTER


A 3-minute reaction period will begin.

Note: A pink color will develop if chlorine is present.

CHLORINE, TOTAL, continued



9. After the timer beeps, place the prepared sample into the cell holder. Tightly cover the sample cell with the instrument cap.



10. Press: **READ**


The cursor will move to the right, then the result in mg/L total chlorine will be displayed.

Note: If the sample temporarily turns yellow after sample addition, or the display flashes "limit", it is due to high chlorine levels. Dilute a fresh sample and repeat the test. A slight loss of chlorine may occur during dilution. Multiply the result by the dilution factor; see Section 1 of the DR/850 Procedures Manual.

Note: Standard Adjust may be performed using a prepared standard (see Standard Adjust in Section 1 of the DR/850 Procedures Manual).

CHLORINE, TOTAL, continued

Using AccuVac Ampuls




1. Enter the stored program number for total chlorine (Cl_2) AccuVac Ampuls.

Press: **PRGM**

The display will show:


PRGM ?

Note: For most accurate results, perform a Reagent Blank Correction using deionized water (see Section 1 of the DR/850 Procedures Manual).




2. Press: **11 ENTER**

The display will show **mg/L, Cl2** and the **ZERO** icon.




3. Fill a sample cell with at least 10 mL of sample (the blank). Collect at least 40 mL of sample in a 50-mL beaker.

Note: Samples must be analyzed immediately and cannot be preserved for later analysis.



4. Place the blank into the cell holder. Tightly cover the sample cell with the instrument cap.




5. Press: **ZERO**

The cursor will move to the right, then the display will show:


0.00 mg/L Cl2

Note: If Reagent Blank Correction is on, the display may flash "limit". See Section 1 of the DR/850 Procedures Manual.




6. Fill a DPD Total Chlorine Reagent AccuVac Ampul with sample.

Note: Keep the tip immersed while the ampul fills completely.



7. Quickly invert the ampul several times to mix. Wipe off any liquid or fingerprints.

Note: A pink color will form if chlorine is present.




8. Press:


TIMER ENTER

A 3-minute reaction period will begin.

CHLORINE, TOTAL, continued



9. When the timer beeps, place the AccuVac Ampul into the cell holder. Tightly cover the ampul with the instrument cap.



10. Press: **READ**

The cursor will move to the right, then the result in mg/L total chlorine will be displayed.

Note: If the sample temporarily turns yellow after sample addition, or the display shows "limit", it is due to high chlorine levels. Dilute a fresh sample and repeat the test. A slight loss of chlorine may occur during dilution. Multiply the result by the appropriate dilution factor; see Section 1 of the DR/850 Procedures Manual.

Note: Standard Adjust may be performed using a prepared standard (see Section 1 of the DR/850 Procedures Manual).

Sampling and Storage

Analyze samples for chlorine **immediately** after collection. Chlorine is a strong oxidizing agent, and it is unstable in natural waters. It reacts rapidly with various inorganic compounds and more slowly oxidizes organic compounds. Many factors, including reactant concentrations, sunlight, pH, temperature, and salinity influence decomposition of chlorine in water.

Avoid plastic containers since these may have a large chlorine demand. **Pretreat glass** sample containers to remove any chlorine demand by soaking in a dilute bleach solution (1 mL commercial bleach to 1 liter of deionized water) for at least 1 hour. Rinse thoroughly with deionized or distilled water. If sample containers are rinsed thoroughly with deionized or distilled water after use, only occasional pre-treatment is necessary.

CHLORINE, TOTAL, continued

Do not use the same sample cells for free and total chlorine. If trace iodide from the total chlorine reagent is carried over into the free chlorine determination, monochloramine will interfere. It is best to use separate, dedicated sample cells for free and total chlorine determinations.

A common error in testing for chlorine is introduced when a representative sample is not obtained. If sampling from a tap, let the water flow for at least 5 minutes to ensure a representative sample. Let the container overflow with the sample several times, then cap the sample containers so there is no headspace (air) above the sample. If sampling with a sample cell, rinse the cell several times with the sample, then carefully fill to the 10-mL mark. Perform the chlorine analysis immediately.

Accuracy Check

Standard Additions Method (using powder pillows)

- a) Snap the top off a LR Chlorine PourRite Ampule Standard Solution.
- b) Use a TenSette Pipet to add 0.1 mL of the standard to the reacted sample (this is the spiked sample). Swirl to mix.
- c) Re-zero the instrument using the original sample (the blank).
- d) Place the spiked sample into the cell holder and press **READ**. Record the results.
- e) Calculate the concentration of mg/L chlorine added to the sample:

$$\text{mg/L chlorine added} = \frac{0.1 (\text{vol. standard added}) \times \text{Certificate value (mg/L Cl}_2\text{)}}{10.1 (\text{sample} + \text{standard volume})}$$

- f) The spiked sample result (*Step d*) should reflect the analyzed sample result + the calculated mg/L Cl₂ added (*Step e*).
- g) If this increase does not occur, see *Standard Additions* in *Section 1* of the *DR/850 Procedures Manual* for more information.

CHLORINE, TOTAL, continued

Standard Additions Method (using AccuVac Ampuls)

- a) Snap the top off a LR Chlorine PourRite Ampule Standard Solution.
- b) Use a graduated cylinder to measure 25 mL of sample into each of two beakers.
- c) Use a TenSette Pipet to add 0.2 mL of the standard to one of the beakers (this is the spiked sample). Swirl to mix.
- d) Fill a DPD Total Chlorine AccuVac completely from each beaker.
- e) Analyze the spiked and unspiked sample as described in the procedure.
- f) Calculate the concentration of mg/L chlorine added to the sample:

$$\text{mg/L chlorine added} = \frac{0.2 \text{ (vol. standard added)} \times \text{Certificate value (mg/L Chlorine)}}{25.2 \text{ (sample + standard volume)}}$$

- g) The spiked sample result should reflect the analyzed sample result + the calculated mg/L Cl₂ added (*Step f*).
- h) If this increase does not occur, see *Standard Additions in Section 1* of the *DR/850 Procedures Manual* for more information.

Method Performance

Precision

In a single laboratory, using a standard solution of 1.00 mg/L chlorine and two lots of reagents with the instrument, a single operator obtained standard deviations of ± 0.01 mg/L chlorine.

In a single laboratory, using a standard solution of 1.00 mg/L chlorine and two representative lots of AccuVac Ampuls with the instrument, a single operator obtained a standard deviation of ± 0.01 mg/L chlorine.

Estimated Detection Limit (EDL)

The estimated detection limit for programs 9 and 11 is 0.02 mg/L Cl₂. For more information on derivation and use of Hach's estimated detection limit, see *Section 1* of the *DR/850 Procedures Manual*.

CHLORINE, TOTAL, continued

Interferences

Interfering Substance	Interference Level and Treatment
Acidity	Greater than 150 mg/L CaCO ₃ . May not develop full color or color may fade instantly. Neutralize to pH 6-7 with 1 N sodium hydroxide. Determine amount to be added on separate sample aliquot, then add the same amount to the sample being tested. Correct for volume addition (see <i>Section 1 of the DR/850 Procedures Manual, Correcting for Volume Additions</i>).
Alkalinity	Greater than 250 mg/L CaCO ₃ . May not develop full color or color may fade instantly. Neutralize to pH 6-7 with 1 N sulfuric acid. Determine amount to be added on separate sample aliquot, then add the same amount to the sample being tested. Correct for volume addition (See <i>Correcting for Volume Additions, Section 1 of the DR/850 Procedures Manual</i>).
Bromine	Interferes at all levels
Chlorine Dioxide	Interferes at all levels
Chloramines, organic	May interfere
Hardness	No effect at less than 1,000 mg/L as CaCO ₃
Iodine	Interferes at all levels
Manganese, Oxidized (Mn ⁴⁺ , Mn ⁷⁺) or Chromium, Oxidized (Cr ⁶⁺)	<ol style="list-style-type: none">1. Adjust sample pH to 6-7.2. Add 3 drops potassium iodide (30 g/L) to a 25-mL sample.3. Mix and wait 1 minute.4. Add 3 drops sodium arsenite (5 g/L) and mix.5. Analyze 10 mL of the treated sample as described in the procedure.6. Subtract the result from this test from the original analysis to obtain the correct chlorine concentration.
Ozone	Interferes at all levels
Peroxides	May interfere
Extreme sample pH and highly buffered samples	Adjust to pH 6-7. See <i>Interferences, Section 1 of the DR/850 Procedures Manual</i> .

CHLORINE, TOTAL, continued

Summary of Method

Chlorine can be present in water as free available chlorine and as combined available chlorine. Both forms can exist in the same water and be determined together as the total available chlorine. Free chlorine is present as hypochlorous acid and/or hypochlorite ion. Combined chlorine exists as monochloramine, dichloramine, nitrogen trichloride, and other chloro derivatives.

The combined chlorine oxidizes iodide in the reagent to iodine. The iodine reacts with DPD (N, N-diethyl-p-phenylenediamine) along with free chlorine present in the sample to form a red color which is proportional to the total chlorine concentration. To determine the concentration of combined chlorine, run free chlorine and total chlorine tests. Subtract the results of the free chlorine test from the results of the total chlorine test to obtain combined chlorine.

Pollution Prevention and Waste Management

Samples treated with sodium arsenite for manganese or chromium interferences will be hazardous wastes as regulated by Federal RCRA for arsenic (D004). See *Section 3* of the *DR/850 Procedures Manual* for more information on proper disposal of these materials.

CHLORINE, TOTAL, continued

REQUIRED REAGENTS & APPARATUS (Using Powder Pillows)

Description	Quantity Required		
	Per Test	Unit	Cat. No.
DPD Total Chlorine Reagent Powder Pillows	1 pillow	100/pkg.....	21056-69
Sample Cell, 10-20-25 mL, with caps	2.....	6/pkg.....	24019-06

REQUIRED REAGENTS & APPARATUS (Using AccuVac Ampuls)

DPD Total Chlorine Reagent AccuVac Ampuls	1 ampul	25/pkg.....	25030-25
Beaker, 50 mL	1.....	each.....	500-41H

OPTIONAL REAGENTS

Chlorine Standard Solution, PourRite Ampule, 25-30 mg/L Cl ₂	20/pkg.....	26300-20
DPD Total Chlorine Reagent, Swiftest Dispenser	250 tests.....	28024-00
Potassium Iodide Solution, 30 g/L.....	100 mL* MDB.....	343-32
Sodium Arsenite, 5 g/L.....	100 mL* MDB.....	1047-32
Sodium Hydroxide Standard Solution, 1.00 N	100 mL* MDB.....	1045-32
Sulfuric Acid Standard Solution, 1.000 N	100 mL* MDB.....	1270-32
Water, deionized.....	4 L.....	272-56

OPTIONAL APPARATUS

AccuVac Snapper	each.....	24052-00
PourRite Ampule Breaker.....	each.....	24846-00
Cylinder, graduated, 25 mL	each.....	508-40
pH Indicator Paper, 1 to 11 pH units	5 rolls/pkg.....	391-33
pH Meter, portable	each.....	HQ11d
pH Electrode, gel, standard, with 1 meter cable	each.....	PHC10101
Pipet, TenSette, 0.1 to 1.0 mL.....	each.....	19700-01
Pipet Tips, for 19700-01 TenSette Pipet	50/pkg.....	21856-96

For Technical Assistance, Price and Ordering

In the U.S.A.—Call 800-227-4224


Outside the U.S.A.—Contact the Hach office or distributor serving you.

* Marked Dropper Bottle - contact Hach for larger sizes.

NITRATE, Low Range (0 to 0.50 mg/L NO₃⁻-N)

Cadmium Reduction Method

For water, wastewater, and seawater*




1. Enter the stored program number for low range nitrate-nitrogen (NO₃⁻-N).

Press: **PRGM**

The display will show:

PRGM ?


Note: For most accurate results, perform a Reagent Blank Correction using deionized water (see Section 1 of the DR/850 Procedures Manual).



2. Press: **55 ENTER**


The display will show **mg/L, NO₃-N** and the **ZERO** icon.

*Note: For alternate forms (NO₃), press the **CONC** key.*



3. Fill a 25-mL graduated mixing cylinder to the 15-mL mark with sample.

Note: Adjust the pH of stored samples before analysis.




4. Add the contents of one NitraVer 6 Nitrate Reagent Powder Pillow to the cylinder. Stopper.

*Note: It is necessary to remove **all** the powder from the foil pillow. Tap the pillow until no more powder pours out. Be sure to remove powder from the corners of the pillow.*

* Seawater requires a manual calibration; see *Interferences*.

NITRATE, Low Range, continued




5. Press:

TIMER ENTER

A 3-minute reaction period will begin. Shake the cylinder vigorously throughout this 3-minute period.

Note: Shaking time and technique influence color development. For most accurate results, analyze a standard solution several times and adjust the shaking time to obtain the correct result.




6. When the timer beeps, the display will show: **2:00 TIMER 2**

Press: **ENTER**


A 2-minute reaction period will begin.

Note: A deposit will remain after the powder dissolves and will not affect results.



7. When the timer beeps, pour 10 mL of the sample into a sample cell.

Note: Do not transfer any cadmium particles.



8. Add the contents of one NitriVer 3 Nitrite Reagent Powder Pillow to the sample cell (the prepared sample). Cap the cell and shake gently for 30 seconds.

Note: A pink color will form if nitrate is present.




9. The display will show: **15:00 TIMER 3**


Press: **ENTER**

A 15-minute reaction period will begin.

Fill another sample cell (the blank) with 10 mL of sample.



10. When the timer beeps, place the blank into the cell holder. Tightly cover the sample cell with the instrument cap.




11. Press: **ZERO**

The cursor will move to the right, then the display will show:


0.00 mg/L NO₃-N

Note: If Reagent Blank Correction is on, the display may flash "limit". See Section 1 of the DR/850 Procedures Manual.



12. Place the prepared sample into the cell holder. Tightly cover the sample cell with the instrument cap.

NITRATE, Low Range, continued



13. Press: **READ**

The cursor will move to the right, then the result in mg/L NO₃⁻-N (or alternate form) will be displayed.

Note: Standard Adjust may be performed using a prepared standard (see Section 1 of the DR/850 Procedures Manual).

Note: Rinse the sample cell and cylinder immediately after use to remove all cadmium particles.

Note: See Pollution Prevention and Waste Management for proper disposal of cadmium.

Sampling and Storage

Collect samples in clean plastic or glass bottles. Store at 4 °C (39 °F) or lower if the sample is to be analyzed within 24 to 48 hours. Warm to room temperature before running the test. For longer storage periods, adjust sample pH to 2 or less with sulfuric acid, ACS (about 2 mL per liter). Sample refrigeration is still required.

Before testing the stored sample, warm to room temperature and neutralize with 5.0 N Sodium Hydroxide Standard Solution. Do not use mercury compounds as preservatives. Correct the test result for volume additions; see *Correction for Volume Additions* (Section 1 of the DR/850 Procedures Manual) for more information.

NITRATE, Low Range, continued

Accuracy Check

Standard additions Method

- a) Fill three 25-mL graduated mixing cylinders with 15 mL of sample.
- b) Snap the neck off a Nitrate-Nitrogen Ampule Standard Solution, 12.0 mg/L NO_3^- -N.
- c) Using the TenSette Pipet, add 0.1, 0.2, and 0.3 mL of the standard to the three samples. Stopper and mix well.
- d) Analyze each sample as described above. The nitrate-nitrogen concentration should increase 0.08 mg/L for each 0.1 mL of standard added.
- e) If these increases do not occur, see *Standard Additions (Section 1 of the DR/850 Procedures Manual)* for more information.

Standard Solution Method

Prepare a 0.20 mg/L nitrate-nitrogen standard by diluting 2.00 mL of a 10.0 mg/L Nitrate-Nitrogen Standard Solution to 100.0 mL with deionized water. Use this standard in place of sample in *Step 3* on page 55.

Standard Adjust

To adjust the calibration curve using the reading obtained with the 0.20-mg/L standard solution, press the **SETUP** key and scroll (using the arrow keys) to the STD setup option. Press **ENTER** to activate the standard adjust option. Then enter **0.20** to edit the standard concentration to match that of the standard used. Press **ENTER** to complete the curve adjustment. If you are using a reagent blank correction, enter the blank correction before the Standard Adjust feature is entered. See *Standard Curve Adjustment, Section 1* of the *DR/850 Procedures Manual* for more information.

NITRATE, Low Range, continued

Method Performance

Precision

In a single laboratory using a standard solution of 0.25 mg/L nitrate-nitrogen (NO_3^- -N) and two representative lots of reagent with the instrument, a single operator obtained a standard deviation of ± 0.03 mg/L nitrate-nitrogen.

Estimated Detection Limit

The estimated detection limit for program 55 is 0.01 mg/L NO_3^- -N. For more information on the estimated detection limit, see *Section 1* of the *DR/850 Procedures Manual*.

Interferences

Interfering Substance	Interference Levels and Treatments
Calcium	100 mg/L
Chloride	Chloride concentrations above 100 mg/L will cause low results. The test may be used at high chloride concentrations (seawater) but a calibration must be done using standards spiked to the same chloride concentration.
Ferric iron	All levels
Nitrite	All levels: This method measures both the nitrate and nitrite in the sample. If nitrite is present, the nitrite-nitrogen test Program 60 should be done on the sample. Pretreat the nitrate-nitrogen sample with the following pretreatment. Then subtract the amount of nitrite found from the results of the LR nitrate-nitrogen test using the pretreated sample. 1. Add 30-g/L Bromine Water dropwise to the sample in <i>Step 3</i> on page 55 until a yellow color remains. Mix after each drop. 2. Add one drop of 30-g/L Phenol Solution to destroy the yellow color. 3. Proceed with the LR Nitrate procedure.
pH	Highly buffered samples or extreme sample pH may exceed the buffering capacity of the reagents and require sample pretreatment.
Strong oxidizing and reducing substances	Interfere at all levels

NITRATE, Low Range, continued

Summary of Method

Cadmium metal reduces nitrates present in the sample to nitrite. The nitrite ion reacts in an acidic medium with sulfanilic acid to form an intermediate diazonium salt which couples to chromotropic acid to form a pink-colored product.

Pollution Prevention and Waste Management

NitaVer 6 Nitrate Reagent contains cadmium metal. Both samples and reagent blanks will contain cadmium (D006) at a concentration regulated as hazardous wastes by the Federal RCRA. Do not pour these solutions down the drain. See *Section 3* of the *DR/850 Procedures Manual* for more information on proper disposal of these materials.

REQUIRED REAGENTS

Low Range Nitrate-Nitrogen Reagent Set 24298-00
Includes: (1) 21071-69, (1) 21072-49

Description	Quantity Required		
	Per Test	Unit	Cat. No.
NitriVer 3 Nitrite Reagent Powder Pillows.....	1 pillow	100/pkg	21071-69
NitraVer 6 Nitrate Reagent Powder Pillows	1 pillow	100/pkg	21072-49

REQUIRED APPARATUS

Cylinder, graduated, mixing, 25 mL 1 each 1896-40
Sample Cell, 10-20-25 mL, with cap 2 6/pkg 24019-06

OPTIONAL REAGENTS

Bromine Water, 30 g/L 29 mL* 2211-20
Nitrate-Nitrogen Standard Solution, 10.0 mg/L as NO_3^- -N 500 mL 307-49
Nitrate-Nitrogen Standard Solution, 12 mg/L as NO_3^- -N,
 Voluette Ampule, 10 mL 16/pkg 14333-10
Phenol Solution, 30 g/L 29 mL 2112-20
Pretreatment Kit, contains: (1) 2112-20, (1) 2211-20 each 2268-00
Sodium Hydroxide Standard Solution, 5.0 N 50 mL* SCDB 2450-26
Sulfuric Acid, ACS 500 mL* 979-49
Water, deionized 4 L 272-56

* Contact Hach for larger sizes

NITRATE, Low Range, continued

OPTIONAL APPARATUS

Description	Unit	Cat. No.
Dropper, for 29-mL bottle	each	2258-00
Flask, volumetric, Class A, 100 mL.....	each	14574-42
pH Indicator Paper, 1 to 11 pH.....	5-roll/pkg	391-33
pH Meter, portable	each	HQ11d
pH Electrode, gel, standard, with 1 meter cable	each	PHC10101
Pipet, serological, 2 mL.....	each	532-36
Pipet, TenSette, 0.1 to 1.0 mL	each	19700-01
Pipet Tips, for 19700-01 TenSette Pipet	50/pkg	21856-96
Pipet, volumetric, Class A, 2.00 mL	each	14515-36
Pipet Filler, safety bulb.....	each	14651-00
PourRite Ampule Breaker	each	24846-00
Thermometer, -10 to 110 °C	each	1877-01

Nitrate at these levels can be determined directly using the Nitrate Ion Selective Electrode (Cat. No. 50235-00).


For Technical Assistance, Price and Ordering

In the U.S.A.—Call 800-227-4224

Outside the U.S.A.—Contact the Hach office or distributor serving you.

NITRITE, Low Range (0 to 0.350 mg/L NO₂⁻-N) For water, wastewater, seawater**Diazotization Method* (Powder Pillows or AccuVac Ampuls);**

USEPA approved for reporting wastewater and drinking water analyses.




1. Enter the stored program number for nitrite-nitrogen (NO₂⁻-N), powder pillows.

Press: **PRGM**

The display will show:

PRGM ?

Note: For most accurate results, perform a Reagent Blank Correction using deionized water (see Section 1 of the DR/850 Procedures Manual).




2. Press: **60 ENTER**

The display will show **mg/L, NO₂-N** and the **ZERO** icon.

*Note: For alternate forms (NO₂⁻, NaNO₂), press the **CONC** key.*



3. Fill a sample cell with 10 mL of sample.




4. Add the contents of one NitriVer 3 Nitrite Reagent Powder Pillow to the sample cell. Cap the cell and shake to dissolve.

Note: Accuracy is not affected by undissolved powder.

* Federal Register, 44(85) 25505 (May 1, 1979)

NITRITE, Low Range, continued



5. Press:


TIMER ENTER

A 15-minute reaction period will begin.

Note: A pink color will develop if nitrite is present.




6. When the timer beeps, fill an empty sample cell with 10 mL of sample (the blank).



7. Wipe the outside of the sample cell with a towel. Place the blank into the cell holder. Tightly cover the sample cell with the instrument cap.

Note: Wiping with a damp cloth, followed by a dry one, removes fingerprints and other marks.




8. Press: **ZERO**


The cursor will move to the right, then the display will show:

0.000 mg/L NO₂-N

Note: If Reagent Blank Correction is on, the display may flash "limit." See Section 1 of the DR/850 Procedures Manual.



9. Place the prepared sample into the cell holder. Tightly cover the sample cell with the instrument cap.




10. Press: **READ**

The cursor will move to the right, then the result in mg/L nitrite-nitrogen (or an alternate form) will be displayed.

NITRITE, Low Range, continued

Using AccuVac Ampuls




1. Enter the stored program number for nitrite-nitrogen (NO_2^- -N), AccuVac Ampuls.

Press: **PRGM**

The display will show:

PRGM ?


Note: For most accurate results, perform a Reagent Blank Correction using deionized water (see Section 1 of the DR/850 Procedures Manual).



2. Press: **62 ENTER**


The display will show **mg/L, NO₂-N** and the **ZERO** icon.

*Note: For alternate forms (NO_2^- , NaNO_2), press the **CONC** key.*




3. Collect at least 40 mL of sample in a 50-mL beaker. Fill a NitriVer 3 Nitrite AccuVac Ampul with the sample.

Note: Keep the tip immersed while the ampul fills completely.



4. Quickly invert the ampul several times to mix. Wipe off any liquid or fingerprints.

Note: Accuracy is not affected by undissolved powder.




5. Press: **TIMER ENTER**

A 15-minute reaction period will begin.


Note: A pink color will develop if nitrite is present.



6. When the timer beeps, fill a sample cell with at least 10 mL of sample (the blank).



7. Place the blank into the cell holder. Tightly cover the sample cell with the instrument cap.




8. Press: **ZERO**

The cursor will move to the right, then the display will show:


0.000 mg/L NO₂-N

Note: If Reagent Blank Correction is on, the display may flash "limit." See Section 1 of the DR/850 Procedures Manual.

NITRITE, Low Range, continued



9. Place the AccuVac Ampul into the cell holder. Tightly cover the ampul with the instrument cap.



10. Press: **READ**

The cursor will move to the right, then the result in mg/L nitrite-nitrogen will be displayed.

Sampling and Storage

Collect samples in clean plastic or glass bottles.

Store at 4 °C (39 °F) or lower and analyze within 48 hours. Warm to room temperature before running the test.

Do not use acid preservatives.

Remove the suspended solids by filtration.

Accuracy Check

Standard Solution Method

Pipet 5.00 mL of a fresh 250 mg/L NO_2^- -N standard into a 250.0 mL volumetric flask. Dilute to the mark with deionized water. This makes a 5.00-mg/L intermediate standard. To prepare a 0.100-mg/L NO_2^- -N standard solution, dilute 10.00 mL of the 5.00-mg/L intermediate standard to 500 mL in a volumetric flask. Prepare this solution immediately before use.

Run the test using the 0.100 mg/L NO_2^- -N standard in place of the sample. Results should be between 0.090 and 0.110 mg/L NO_2^- -N.

NITRITE, Low Range, continued

Method Performance

Precision

In a single laboratory, using a standard solution of 0.250 mg/L nitrite-nitrogen and two representative lots of reagent with the instrument, a single operator obtained a standard deviation of ± 0.001 mg/L NO_2^- -N for the powder pillow method and ± 0.003 mg/L NO_2^- -N for the AccuVac method.

Estimated Detection Limit

The estimated detection limit for programs 60 and 62 is 0.005 mg/L NO_2^- -N. For more information on derivation and use of Hach's estimated detection limit, see *Section 1* of the *DR/850 Procedures Manual*.

Interferences

Interfering Substance	Interference Levels
Antimonious ions	Interfere by causing precipitation
Auric ions	Interfere by causing precipitation
Bismuth ions	Interfere by causing precipitation
Chloroplatinate ions	Interfere by causing precipitation
Cupric ions	Cause low results
Ferric ions	Interfere by causing precipitation
Ferrous ions	Cause low results
Lead ions	Interfere by causing precipitation
Mercurous ions	Interfere by causing precipitation
Metavanadate ions	Interfere by causing precipitation
Nitrate	Very high levels of nitrate (> 100 mg/L nitrate as N) appear to undergo a slight amount of reduction to nitrite, either spontaneously or during the course of the test. A small amount of nitrite will be found at these levels.
Silver ions	Interfere by causing precipitation
Strong oxidizing and reducing substances	Interfere at all levels

Summary of Method

Nitrite in the sample reacts with sulfanilic acid to form an intermediate diazonium salt. This couples with chromotropic acid to produce a pink colored complex directly proportional to the amount of nitrite present.

NITRITE, Low Range, continued

REQUIRED REAGENTS

Description	Quantity Required		Unit	Cat. No.
	Per Test			
NitriVer 3 Nitrite Reagent Powder Pillows.....	1 pillow	100/pkg.....	21071-69
or				
NitriVer 3 Nitrite Reagent AccuVac Ampuls.....	1 ampul	25/pkg.....	25120-25

REQUIRED APPARATUS

Beaker, 50 mL (for AccuVac procedure).....	1	each.....	500-41H
or				
Sample Cells, 10-20-25 mL (powder pillow procedure)2			6/pkg.....	24019-06

OPTIONAL REAGENTS

Nitrite Standard Solution, 250 mg/L			500 mL	23402-49
Water, deionized.....			4 L.....	272-56

OPTIONAL APPARATUS

AccuVac Snapper			each.....	24052-00
Flask, volumetric, 250 mL			each.....	14574-46
Flask, volumetric, 500 mL			each.....	14574-49
Pipet, serological, 10 mL			each.....	532-38
Pipet, TenSette, 1 to 10 mL.....			each.....	19700-01
Pipet Tips for 19700-01 TenSette Pipet			50/pkg.....	21856-96
Pipet, volumetric, Class A, 5.00 mL.....			each.....	14515-37
Pipet, volumetric, Class A, 10.00 mL.....			each.....	14515-38
Pipet Filler, safety bulb			each.....	14651-00
Thermometer, -10 to 110 °C.....			each.....	1877-01

For Technical Assistance, Price and Ordering


In the U.S.A.—Call 800-227-4224

Outside the U.S.A.—Contact the Hach office or distributor serving you.

NITROGEN, AMMONIA (0 to 0.50 mg/L NH₃-N)

Salicylate Method*

For water, wastewater, and seawater




1. Enter the stored program number for ammonia-nitrogen (NH₃-N).

Press: **PRGM**

The display will show:


PRGM ?




2. Press: **64 ENTER**

The display will show **mg/L, NH₃-N** and the **ZERO** icon.


*Note: For alternate forms (NH₃, NH₄), press the **CONC** key.*




3. Fill a sample cell with 10 mL of deionized water (the blank).



4. Fill a second sample cell with 10 mL of the sample.




5. Add the contents of one Ammonia Salicylate Reagent Powder Pillow to each sample cell. Cap both cells and shake to dissolve.




6. Press: **TIMER ENTER**

A 3-minute reaction period will begin.



7. After the timer beeps add the contents of one Ammonia Cyanurate Reagent Powder Pillow to each sample cell. Cap the cells and shake to dissolve the reagent.

Note: A green color will develop if ammonia-nitrogen is present.




8. The display will show: **15:00 TIMER 2**

Press: **ENTER**


A 15-minute reaction period will begin.

* Adapted from *Clin. Chim. Acta.*, 14 403 (1966).


NITROGEN, AMMONIA, continued




9. After the timer beeps, place the blank into the cell holder. Tightly cover the sample cell with the instrument cap.



10. Press: **ZERO**
The cursor will move to the right, then the display will show:
0.00 mg/L NH₃-N



11. Place the prepared sample into the cell holder. Tightly cover the sample cell with the instrument cap.



12. Press: **READ**
The cursor will move to the right, then the result in mg/L ammonia-nitrogen will be displayed.

Note: Standard Adjust may be performed using a prepared standard (see Standard Adjust in Section 1 of the DR/850 Procedures Manual).

Sampling and Storage

Collect samples in clean plastic or glass bottles. Most reliable results are obtained when samples are analyzed as soon as possible after collection.

If chlorine is known to be present, the sample must be treated immediately with sodium thiosulfate. Add one drop of Sodium Thiosulfate Standard Solution, 0.1 N, for each 0.3 mg of chlorine present in a 1-liter sample.

To preserve the sample, adjust the pH to 2 or less with concentrated sulfuric acid (about 2 mL per liter). Store samples at 4 °C or less. Samples preserved in this manner can be stored up to 28 days. Just before testing the stored sample, warm to room temperature and neutralize with 5.0 N Sodium Hydroxide Standard Solution. Correct the test result for volume additions; see *Correction for Volume Additions*, in Section 1 of the *DR/850 Procedures Manual* for more detailed information.

Accuracy Check

Standard Additions Method

- a) Fill three 25-mL mixing cylinders with 20 mL of sample.
- b) Use the TenSette Pipet to add 0.1, 0.2, and 0.3 mL of Ammonium Nitrogen Standard, 10 mg/L as $\text{NH}_3\text{-N}$ to the three samples. Stopper the cylinders and mix well.
- c) Analyze a 10-mL portion of sample as described above. The ammonia-nitrogen concentration should increase 0.05 mg/L for each 0.1 mL of standard added.
- d) If these increases do not occur, see *Standard Additions (Section 1 of the DR/850 Procedures Manual)* for more information.

Standard Solution Method

Prepare a 0.40 mg/L ammonia-nitrogen standard by diluting 4.00 mL of the Ammonia-Nitrogen Standard Solution, 10 mg/L, to 100 mL with deionized water. Or, using the TenSette Pipet, prepare a 0.40 mg/L ammonia-nitrogen standard by diluting 0.8 mL of a Ammonia-Nitrogen Voluette Standard Solution, 50 mg/L as $\text{NH}_3\text{-N}$, to 100 mL with deionized water.

Method Performance

Precision

In a single laboratory using a standard solution of 0.40 mg/L ammonia-nitrogen ($\text{NH}_3\text{-N}$) and two representative lots of reagent with the instrument, a single operator obtained a standard deviation of ± 0.02 mg/L ammonia-nitrogen.

Estimated Detection Limit

The estimated detection limit for program 64 is 0.02 mg/L $\text{NH}_3\text{-N}$. For more information on the estimated detection limit, see *Section 1 of the DR/850 Procedures Manual*.

NITROGEN, AMMONIA, continued

Interferences.

Interfering Substance	Interference Level and Treatments
Calcium	Greater than 1000 mg/L as CaCO ₃
Glycine, hydrazine	Less common. Will cause intensified colors in the prepared sample.
Iron	All levels. Correct for iron interference as follows: 1. Determine the amount of iron present in the sample using one of the Total Iron procedures. 2. Prepare a deionized water sample containing the same iron concentration as the original sample. Run the procedure on this solution to determine the interference due to iron. Subtract this value from the result in <i>Step 12</i> on page 70 obtained on the original sample.
Magnesium	Greater than 6000 mg/L as CaCO ₃
Nitrate	Greater than 100 mg/L as NO ₃ ⁻ -N
Nitrite	Greater than 12 mg/L as NO ₂ ⁻ -N
Phosphate	Greater than 100 mg/L as PO ₄ ³⁻ -P
Sulfate	Greater than 300 mg/L as SO ₄ ²⁻
Sulfide	Sulfide will intensify the color. Eliminate sulfide interference as follows: 1. Measure about 350 mL of sample in a 500-mL Erlenmeyer flask. 2. Add the contents of one Sulfide Inhibitor Reagent Powder Pillow. Swirl to mix. 3. Filter the sample through a folded filter paper. 4. Use the filtered solution in <i>Step 3</i> on page 69.
Turbidity, sample color	Turbidity and sample color will give erroneous high values. Samples with severe interferences require distillation. Albuminoid nitrogen samples also require distillation. Hach recommends the distillation procedure using the Hach General Purpose Distillation Set. See optional apparatus listing.

Summary of Method

Ammonia compounds combine with chlorine to form monochloramine. Monochloramine reacts with salicylate to form 5-aminosalicylate. The 5-aminosalicylate is oxidized in the presence of a sodium nitroprusside catalyst to form a blue-colored compound. The blue color is masked by the yellow color from the excess reagent present to give a final green-colored solution.

NITROGEN, AMMONIA, continued

REQUIRED REAGENTS AND APPARATUS

Ammonia-Nitrogen Reagent Set for 10-mL samples (100 Tests).....			Cat. No.
			26680-00
Includes: (2) 26531-99, (2) 26532-99			

Description	Quantity Required		Unit	Cat. No.
	Per Test			
Ammonia Cyanurate Reagent Powder Pillows	2 pillows	100/pkg	26531-99	
Ammonia Salicylate Reagent Powder Pillows	2 pillows	100/pkg	26532-99	
Sample Cell, 10-20-25 mL, with cap	2	6/pkg	24019-06	

OPTIONAL REAGENTS

Ammonia-Nitrogen Standard Solution, 10 mg/L as NH ₃ -N	500 mL	153-49
Ammonia-Nitrogen, PourRite Ampule, 50 mg/L as NH ₃ -N, 2 mL	20/pkg	14791-20
Cylinder, graduated, mixing, 25 mL	each	20886-40
Sodium Hydroxide Standard Solution, 1.00 N.....	100 mL MDB	1045-32
Sodium Hydroxide Standard Solution, 5.0 N.....	50 mL SCDB	2450-26
Sodium Thiosulfate Standard Solution, 0.1 N.....	100 mL MDB	323-32
Sulfide Inhibitor Reagent Powder Pillows	100/pkg	2418-99
Sulfuric Acid, concentrated, ACS	500 mL	979-49
Sulfuric Acid Standard Solution, 1.000 N.....	100 mL MDB	1270-32
Water, deionized	4 L	272-56

OPTIONAL APPARATUS

Cylinder, graduated, polypropylene, 500 mL.....	each	1081-49
Distillation Heater and Support Apparatus, 115 Vac	each	22744-00
Distillation Heater and Support Apparatus, 230 Vac	each	22744-02
Distillation Set, General Purpose	each	22653-00
Filter Paper, folded, 12.5 cm	100/pkg	1894-57
Flask, Erlenmeyer, polypropylene, 500 mL	each	1082-49
Flask, volumetric, Class A, 100 mL.....	each	14574-42
Funnel, poly, 65 mm.....	each	1083-67
pH Meter, portable	each	HQ11d
pH Electrode, gel, standard, with 1 meter cable	each	PHC10101
Pipet Filler, safety bulb.....	each	14651-00
Pipet, TenSette, 0.1 to 1.0 mL.....	each	19700-01
Pipet Tips, for 19700-01 TenSette Pipet	50/pkg	21856-96
Pipet, volumetric, Class A, 2.0 mL	each	14515-36
PourRite Ampule Breaker	each	24846-00
Thermometer, -10 to 110 °C.....	each	1877-01

For Technical Assistance, Price and Ordering


In the U.S.A.—Call 800-227-4224

Outside the U.S.A.—Contact the Hach office or distributor serving you.

pH (6.5 to 8.5 pH units)

Colorimetric pH Determination Using Phenol Red

For water and wastewater




1. Enter the stored program number for the pH method.

Press: **PRGM**


The display will show:

PRGM ?




2. Press: **75 ENTER**


The display will show **PH** and the **ZERO** icon.



3. Fill a sample cell with 10 mL of sample (the blank).



4. Place the blank in the cell holder. Tightly cover the sample cell with the instrument cap.



5. Press: **ZERO**


The cursor will move to the right, then the display will show:

6.0 PH




6. Fill another cell with 10 mL of sample.


Note: Sample temperature must be 21 to 29 °C.



7. Using a disposable pipet, add 1 mL of Phenol Red Indicator Solution to the cell (the prepared sample). Cap the sample cell and invert twice to mix.



8. Place the prepared sample into the cell holder. Tightly cover the sample cell with the instrument cap.



9. Press: **READ**

The cursor will move to the right, then the result in pH units will be displayed.

Note: Use of the Standard Adjust feature is highly recommended. See Accuracy Check.

Note: Any reading below 6.5 pH units will be erroneous.

Sampling and Storage

Analyze samples immediately for best results.

Accuracy Check

Standard Solution Method

Using a clear pH 7.0 buffer solution as the sample, perform the pH procedure as described above.

Method Performance

Precision

In a single laboratory using a standard solution of pH 7.0 and two lots of reagent with the instrument, a single operator obtained a standard deviation of less than 0.1 pH units.

Estimated Detection Limit

The estimated detection limit for program 75 is a pH of 6.5.

Standard Adjust

To adjust the calibration curve using the reading obtained with the 7.0 buffer solution, press the **SETUP** key and scroll (using the arrow keys) to the STD setup option. Press **ENTER** to activate the standard adjust option. Then enter **7.0** to edit the standard concentration to match that of the standard used. See *Section 1, Standard Curve Adjustment* of the *DR/850 Procedures Manual* for more information. Press **ENTER** to complete the curve adjustment.

Interferences

Chlorine does not interfere at levels of 6 mg/L or lower.

Salt water (sea water) will interfere and cannot be analyzed using this method.

Summary of Method

This method uses a sulfonphthalein indicator (Phenol Red) to determine pH colorimetrically. Phenol Red has a working range of pH 6.8 (yellow) to 8.2 (red).

REQUIRED REAGENTS & APPARATUS

Description	Quantity Required		Units	Cat. No.
	Per Test			
Dropper, 0.5 and 1.0 mL marks.....	1		20/pkg	21247-20
Phenol Red Indicator Solution, spec grade	1.0 mL		50 mL	26575-12
Sample Cells, 10-20-25 mL, with cap.....	2		6/pkg	24019-06

OPTIONAL REAGENTS

pH 7.0 Buffer Solution			500 mL	12222-49
------------------------------	--	--	--------	----------


OPTIONAL APPARATUS

Thermometer, -10 to 110 °C.....			each	1877-01
---------------------------------	--	--	------	---------

For Technical Assistance, Price and Ordering

In the U.S.A.—Call 800-227-4224

Outside the U.S.A.—Contact the Hach office or distributor serving you.

PHOSPHORUS, REACTIVE (0 to 2.50 mg/L PO₄³⁻) For water, wastewater, seawater(Also called Orthophosphate) **PhosVer 3 (Ascorbic Acid) Method*****(Powder Pillows or AccuVac Ampuls)****USEPA Accepted for wastewater analysis reporting******Using Powder Pillows**


1. Enter the stored program number for reactive phosphorus, ascorbic acid method.

Press: **PRGM**

The display will show:

PRGM ?


Note: For most accurate results, perform a Reagent Blank Correction using deionized water (see Section 1 of the DR/850 Procedures Manual).



2. Press: **79 ENTER**

The display will show **mg/L, PO₄** and the **ZERO** icon.


*Note: For alternate forms (P, P₂O₅), press the **CONC** key.*



3. Fill a sample cell with 10 mL of sample.

Note: For samples with extreme pH, see Interferences following these steps.

Note: Clean glassware with 1:1 HCl. Rinse again with deionized water. Do not use detergents containing phosphates to clean glassware.




4. Add the contents of one PhosVer 3 Phosphate Powder Pillow for 10-mL sample to the cell (the prepared sample). Shake for 15 seconds.

Note: A blue color will form if phosphate is present.

* Adapted from *Standard Methods for the Examination of Water and Wastewater*.

** Procedure is equivalent to USEPA method 365.2 and Standard Method 4500-PE for wastewater.

PHOSPHORUS, REACTIVE, continued



5. Press:


TIMER ENTER

A 2-minute reaction period will begin. Perform *Steps 6* through *8* during the reaction period.

Note: If the acid-persulfate digestion was used, an 8 to 10 minute reaction period is required.




6. Fill another sample cell with 10 mL of sample (the blank).



7. Place the blank into the cell holder. Tightly cover the sample cell with the instrument cap.

Press: **EXIT**




8. Press: **ZERO**


The cursor will move to the right, then the display will show:

0.00 mg/L PO₄

Note: If Reagent Blank Correction is on, the display may flash "limit". See Section 1 of the DR/850 Procedures Manual.



9. After the timer beeps, place the prepared sample into the cell holder. Tightly cover the sample cell with the instrument cap.




10. Press: **READ**

The cursor will move to the right, then the result in mg/L phosphate (PO₄³⁻) will be displayed.

Note: Standard Adjust may be performed using a 2.0-mg/L PO₄³⁻-standard; see Section 1 of the DR/850 Procedures Manual.

Using AccuVac Ampuls




1. Enter the stored program number for reactive phosphorus-ascorbic acid method.

Press: **PRGM**

The display will show:

PRGM ?


Note: For most accurate results, perform a Reagent Blank Correction using deionized water (see Section 1 of the DR/850 Procedures Manual).



2. Press: **79 ENTER**

The display will show **mg/L, PO₄** and the **ZERO** icon.


*Note: For alternate forms (P, P₂O₅), press the **CONC** key.*



3. Fill a sample cell (the blank) with at least 10 mL of sample. Collect at least 40-mL of sample in a 50-mL beaker.

Note: For samples with extreme pH, see Interferences.


Note: Clean glassware with 1:1 HCl. Rinse again with deionized water. Do not use detergent containing phosphates to clean glassware.



4. Fill a PhosVer 3 Phosphate AccuVac Ampul with sample.

Note: Keep the tip immersed while the ampul fills completely.


PHOSPHORUS, REACTIVE, continued



5. Place an ampul cap securely over the tip of the ampul. Shake the ampul for about 30 seconds. Wipe off any liquid or fingerprints.

Note: A blue color will form if phosphate is present.

Note: Accuracy is not affected by undissolved powder.




6. Press:

TIMER ENTER


A 2-minute reaction period will begin. Perform Steps 7 and 8 during the reaction period.

Note: Use an 8 to 10 minute reaction period if determining total phosphorus following the acid-persulfate digestion.



7. Place the blank into the cell holder. Tightly cover the sample cell with the instrument cap.

Press: **EXIT**




8. Press: **ZERO**


The cursor will move to the right, then the display will show:

0.00 mg/L PO4

Note: If Reagent Blank Correction is on, the display may flash "limit". See Section 1 of the DR/850 Procedures Manual.



9. After the timer beeps, place the AccuVac Ampul into the cell holder. Tightly cover the ampul with the instrument cap.



10. Press: **READ**

The cursor will move to the right, then the result in mg/L phosphate (PO_4^{3-}) will be displayed.

Note: Standard Adjust may be performed using a 2.0-mg/L PO_4^{3-} standard; see Section 1 of the DR/850 Procedures Manual.

Sampling and Storage

Collect sample in plastic or glass bottles that have been cleaned with 1:1 Hydrochloric Acid Solution and rinsed with deionized water. Do not use commercial detergents containing phosphate for cleaning glassware used in this test.

Analyze samples immediately after collection for best results. If prompt analysis is impossible, preserve samples for up to 48 hours by filtering immediately and storing samples at 4 °C. Warm to room temperature before analysis.

Accuracy Check

Standard Additions Method

- a) Fill three 25-mL graduated mixing cylinders with 25 mL of sample.
- b) Snap the neck off a Phosphate PourRite Ampule Standard Solution, 50 mg/L as PO_4^{3-} .
- c) Use the TenSette Pipet to add 0.1 mL, 0.2 mL, and 0.3 mL of standard, respectively, to the three mixing cylinders. Stopper each and mix thoroughly.
- d) For analysis with AccuVacs, transfer solutions to dry, clean 50-mL beakers to fill the AccuVac Ampuls. For analysis with powder pillows, transfer only 10 mL of solution to the sample cells.
- e) Analyze each standard addition sample as described in the procedure. The phosphate concentration should increase 0.2 mg/L PO_4^{3-} for each 0.1 mL of standard added.
- f) If these increases do not occur, see *Standard Additions in Section 1 of the DR/850 Procedures Manual* for more information.

Standard Solution Method

Prepare a 2.0 mg/L PO_4^{3-} standard solution by pipetting 4.0 mL of Phosphate Standard Solution, 50 mg/L as PO_4^{3-} , into an acid-washed Class A 100-mL volumetric flask. Dilute to volume with deionized water. Stopper and invert to mix. Use this solution in place of the sample in the procedure to insure the accuracy of the test. The mg/L PO_4^{3-} reading should be 2.00 mg/L.

PHOSPHORUS, REACTIVE, continued

Method Performance

Precision

In a single laboratory using a standard solution of 1.00 mg/L PO₄³⁻ and two lots of reagents with the instrument, a single operator obtained a standard deviation of ± 0.05 mg/L PO₄³⁻.

In a single laboratory using a standard solution of 1.00 mg/L PO₄³⁻ and two representative lots of AccuVac Ampuls with the instrument, a single operator obtained a standard deviation of ± 0.03 mg/L PO₄³⁻.

Estimated Detection Limit (EDL)

The EDL for program 79 is 0.05 mg/L PO₄. For more information on the estimated detection limit, see *Section 1* of the *DR/850 Procedures Manual*.

Interference

Interfering Substance	Interference Levels and Treatments
Aluminum	Greater than 200 mg/L
Arsenate	All levels
Chromium	Greater than 100 mg/L
Copper	Greater than 10 mg/L
Hydrogen sulfide	All levels
Iron	Greater than 100 mg/L
Nickel	Greater than 300 mg/L
Silica	Greater than 50 mg/L
Silicate	Greater than 10 mg/L
Turbidity or color	Large amounts may cause inconsistent results in the test because the acid present in the powder pillows may dissolve some of the suspended particles and because of variable desorption of orthophosphate from the particles. For highly turbid or colored samples, add the contents of one Phosphate Pretreatment Pillow to 25 mL of sample. Mix well. Use this solution to zero the instrument.
Zinc	Greater than 80 mg/L
Highly buffered samples or extreme sample pH	May exceed the buffering capacity of the reagents and require sample pretreatment. pH 2 to 10 is recommended.

PHOSPHORUS, REACTIVE, continued

Summary of Method

Orthophosphate reacts with molybdate in an acid medium to produce a Phosphomolybdate complex. Ascorbic acid then reduces the complex, giving an intense molybdenum blue color.

REQUIRED REAGENTS & APPARATUS (Using Powder Pillows)

Description	Quantity Required		Unit	Cat. No.
	Per Test			
PhosVer 3 Phosphate Reagent Powder Pillows, 10 mL sample size	1 Pillow	100/pkg	21060-69	
Sample Cell, 10-20-25 mL, with cap	2	6/pkg	24019-06	

REQUIRED REAGENTS & APPARATUS (Using AccuVac Ampuls)

PhosVer 3 Phosphate Reagent AccuVac Ampuls....	1 ampul	25/pkg	25080-25
Beaker, 50 mL	1	each	500-41H
Cap, ampul, blue	1	25/pkg	1731-25
Sample Cell, 10-20-25 mL, with cap	1	6/pkg	24019-06

OPTIONAL REAGENTS

Drinking Water Mixed Parameter Standard, Inorganic	500 mL	28330-49
Hydrochloric Acid Standard Solution, 6.0 N (1:1)	500 mL	884-49
Phosphate Pretreatment Powder Pillows	50/pkg	14501-99
Phosphate Standard Solution, 1 mg/L	500 mL	2569-49
Phosphate Standard Solution, PourRite Ampule, 50 mg/L, 2 mL	20/pkg	171-20H
Phosphate Standard Solution, Voluette Ampule, 50 mg/L, 10 mL	16/pkg	171-10
Sodium Hydroxide Standard Solution, 5.0 N	100 mL* MDB	2450-32
Wastewater Effluent Mixed Parameter Standard, Inorganic	500 mL	28332-49
Water, deionized	4 L	272-56

* Larger sizes available.

PHOSPHORUS, REACTIVE, continued

OPTIONAL APPARATUS

Description	Unit	Cat. No.
AccuVac Snapper	each.....	24052-00
Ampule Breaker Kit for 10-mL ampules	each.....	21968-00
Aspirator, vacuum pump	each.....	2131-00
Cylinder, graduated, mixing, 25 mL, tall (3 required)	each.....	20886-40
Filter Holder, 47 mm, 300 mL, graduated	each.....	13529-00
Filter Membrane, 47 mm, 0.45 micron	100/pkg.....	13530-00
Flask, filtering, 500 mL	each.....	546-00
Flask, volumetric, Class A, 100 mL	each.....	14574-42
pH Indicator Paper, 1 to 11 pH	5 rolls/pkg	391-33
pH Meter, portable	each.....	HQ11d
pH Electrode, gel, standard, with 1 meter cable	each.....	PHC10101
Pipet, 2 mL, serological	each.....	532-36
Pipet, TenSette, 0.1 to 1.0 mL TenSette Pipet	each.....	19700-01
Pipet Tips, for 19700-01	50/pkg.....	21856-96
Pipet Filler, safety bulb	each.....	14651-00
Pipet, volumetric, Class A, 4.00 mL	each.....	14515-04
PourRite Ampule Breaker	each.....	24846-00


For Technical Assistance, Price and Ordering

In the U.S.A.—Call 800-227-4224

Outside the U.S.A.—Contact the Hach office or distributor serving you.

SULFIDE (0 to 0.70 mg/L S²⁻)

For water, wastewater and seawater


Methylene Blue Method* USEPA accepted for reporting wastewater analysis**

1. Enter the stored program number for sulfide (S).

Press: **PRGM**


The display will show:

PRGM ?



2. Press: **93 ENTER**


The display will show **mg/L, S** and the **ZERO** icon.



3. Pipet 25 mL of sample into a clean sample cell. This will be the prepared sample.

Note: Samples must be analyzed immediately and cannot be preserved for later analysis. Use a pipet to avoid agitation.

Note: For field testing, a 25-mL graduated cylinder may be used.




4. Fill a second sample cell with 25 mL of deionized water (the blank).

* Adapted from *Standard Methods for the Examination of Water and Wastewater*.


** Procedure is equivalent to USEPA method 376.2 or Standard Method 4500-S²⁻ D for wastewater.

SULFIDE, continued




5. Add 1.0 mL of Sulfide 1 Reagent to each cell. Swirl to mix.

Note: Use the calibrated 1-mL dropper.



6. Add 1.0 mL of Sulfide 2 Reagent to each cell. Immediately swirl to mix.


Note: A pink color will develop, then the solution will turn blue if sulfide is present.




7. Press:

TIMER ENTER

A 5-minute reaction period will begin.




8. Place the blank into the cell holder. Tightly cover the sample cell with the instrument cap.




9. Press: **ZERO**

The cursor will move to the right, then the display will show:

0.00 mg/L S



10. After the timer beeps, place the prepared sample into the cell holder. Tightly cover the sample cell with the instrument cap.



11. Press: **READ**

The cursor will move to the right, then the result in mg/L sulfide will be displayed.

Note: Some sulfide loss may occur if dilution is necessary.

Note: Standard Adjust may be performed using a prepared standard (see Standard Adjust in Section 1 of the DR/850 Procedures Manual).

SULFIDE, continued

Sampling and Storage

Collect samples in clean plastic or glass bottles. Fill completely and cap tightly. Avoid excessive agitation or prolonged exposure to air. Analyze samples immediately.

Method Performance

Precision

In a single laboratory, using standard solutions of 0.73 mg/L sulfide and two representative lots of reagent with the instrument, a single operator obtained a standard deviation of ± 0.02 mg/L sulfide.

Estimated Detection Limit (EDL)

The EDL for program 93 is 0.01 mg/L S^{2-} . For more information on derivation and use of Hach's estimated detection limit, see *Section 1* of the *DR/850 Procedures Manual*.

Interferences

Interfering Substance	Interference Levels and Treatments
Strong reducing substances (sulfite, thiosulfate and hydrosulfite)	Interfere by reducing the blue color or preventing its development.
Sulfide, high levels	High concentrations of sulfide may inhibit full color development and require sample dilution. Some sulfide loss may occur when the sample is diluted.
Turbidity	For turbid samples, prepare a sulfide-free blank as follows. Use it in place of the deionized water blank in the procedure. 1. Measure 25 mL of sample into a 50-mL Erlenmeyer flask. 2. Add Bromine Water dropwise with constant swirling until a permanent yellow color just appears. 3. Add Phenol Solution dropwise until the yellow color just disappears. Use this solution in <i>Step 4</i> on page <i>87</i> in place of deionized water.

Soluble Sulfides

Determine soluble sulfides by centrifuging the sample in completely filled, capped tubes and analyzing the supernatant. Insoluble sulfides are then estimated by subtracting the soluble sulfide concentration from the total sulfide result.

SULFIDE, continued

Summary of Method

Hydrogen sulfide and acid-soluble metal sulfides react with N, N-dimethyl-p-phenylenediamine oxalate to form methylene blue. The intensity of the blue color is proportional to the sulfide concentration. High sulfide levels in oil field waters may be determined after dilution.

Pollution Prevention and Waste Management

Sulfide 2 Reagent contains potassium dichromate. The final solution will contain hexavalent chromium (D007) at a concentration regulated as a hazardous waste by Federal RCRA. See *Section 3* of the *DR/850 Procedures Manual* for more information on proper disposal of these materials.

REQUIRED REAGENTS

	Cat. No.
Sulfide Reagent Set (100 tests).....	22445-00
Includes: (2) 1816-32, (2) 1817-32	

Description	Quantity Required		Units	Cat. No.
	Per Test			
Sulfide 1 Reagent.....	2 mL.....	100 mL	MDB.....	1816-32
Sulfide 2 Reagent.....	2 mL.....	100 mL	MDB.....	1817-32
Water, deionized.....	25 mL.....	4 L	272-56

REQUIRED APPARATUS

Cylinder, graduated, 25 mL	1	each.....	508-40
Pipet, volumetric, Class A, 25.00 mL.....	1	each.....	14515-40
Pipet Filler, safety bulb	1	each.....	14651-00
Sample Cell, 10-20-25 mL, with cap.....	2	6/pkg.....	24019-06

OPTIONAL REAGENTS

Bromine Water, 30 g/L.....	29 mL.....	2211-20
Phenol Solution, 30 g/L	29 mL.....	2112-20
Sodium Sulfide, hydrate	114 g.....	785-14

OPTIONAL APPARATUS

Bottle, Wash, 250 mL	each.....	620-31
Dropper, for 29 mL bottle	each.....	2258-00
Flask, Erlenmeyer, 50 mL.....	each.....	505-41
<i>Standard Methods for the Examination of Water and Wastewater</i>	each.....	22708-00

For Technical Assistance, Price and Ordering

In the U.S.A.—Call 800-227-4224

Outside the U.S.A.—Contact the Hach office or distributor serving you.



GENERAL INFORMATION

At Hach Company, customer service is an important part of every product we make.

With that in mind, we have compiled the following information for your convenience.

HOW TO ORDER

U.S.A. Customers

By Telephone:

6:30 a.m. to 5:00 p.m. MST

Monday through Friday

800-227-HACH (800-227-4224)

By Fax:

970-669-2932

By Mail:

Hach Company

P. O. Box 389

Loveland, Colorado 80539-0389 U.S.A.

Ordering information by e-mail: orders@hach.com

Information Required:

- Hach account number (if available)
- Your name and phone number
- Purchase order number
- Brief description or model number
- Billing address
- Shipping address
- Catalog number
- Quantity

International Customers

Hach maintains a worldwide network of dealers and distributors. To locate the representative nearest you, send an e-mail to: intl@hach.com or contact:

Hach Company World Headquarters: Loveland, Colorado, U.S.A.

Telephone: (970) 669-3050; Fax: (970) 669-2932

Technical and Customer Service (U.S.A. only)

Hach Technical and Customer Service Department personnel are eager to answer questions about our products and their use. Specialists in analytical methods, they are happy to put their talents to work for you.

Call 1-800-227-4224 or e-mail techhelp@hach.com

REPAIR SERVICE

Authorization must be obtained from Hach Company before sending any items for repair. Please contact the Hach Service Center serving your location.

In the United States:

Hach Company
Ames Service
100 Dayton Avenue
Ames, Iowa 50010
(800) 227-4224 (USA only)
Telephone: (515) 232-2533
FAX: (515) 232-3835

In Canada:

Hach Sales & Service Canada Ltd.
1313 Border Street, Unit 34
Winnipeg, Manitoba
R3H 0X4
(800) 665-7635 (Canada only)
Telephone: (204) 632-5598
FAX: (204) 694-5134
E-mail: canada@hach.com

In Latin America, the Caribbean, the Far East, Indian Subcontinent, Africa, Europe, or the Middle East:

Hach Company World Headquarters,
P.O. Box 389
Loveland, Colorado 80539-0389 U.S.A.
Telephone (970) 669-3050
FAX: (970) 669-2932
E-mail: intl@hach.com

LIMITED WARRANTY

Hach Company warrants its products to the original purchaser against any defects that are due to faulty material or workmanship for a period of one year from date of shipment unless otherwise noted in the product manual.

In the event that a defect is discovered during the warranty period, Hach Company agrees that, at its option, it will repair or replace the defective product or refund the purchase price excluding original shipping and handling charges. Any product repaired or replaced under this warranty will be warranted only for the remainder of the original product warranty period.

This warranty does not apply to consumable products such as chemical reagents; or consumable components of a product, such as, but not limited to, lamps and tubing.

Contact Hach Company or your distributor to initiate warranty support. Products may not be returned without authorization from Hach Company.

Limitations

This warranty does not cover:

- Damage caused by acts of God, natural disaster, labor unrest, acts of war (declared or undeclared), terrorism, civil strife or acts of any governmental jurisdiction
- Damage caused by misuse, neglect, accident or improper application or installation
- Damage caused by any repair or attempted repair not authorized by Hach Company
- Any product not used in accordance with the instructions furnished by Hach Company
- Freight charges to return merchandise to Hach Company
- Freight charges on expedited or express shipment of warranted parts or product
- Travel fees associated with on-site warranty repair

This warranty contains the sole express warranty made by Hach Company in connection with its products. All implied warranties, including without limitation, the warranties of merchantability and fitness for a particular purpose, are expressly disclaimed.

LIMITED WARRANTY, continued

Some states within the United States do not allow the disclaimer of implied warranties and if this is true in your state the above limitation may not apply to you. This warranty gives you specific rights, and you may also have other rights that vary from state to state.

This warranty constitutes the final, complete, and exclusive statement of warranty terms and no person is authorized to make any other warranties or representations on behalf of Hach Company.

Limitation of Remedies

The remedies of repair, replacement or refund of purchase price as stated above are the exclusive remedies for the breach of this warranty. On the basis of strict liability or under any other legal theory, in no event shall Hach Company be liable for any incidental or consequential damages of any kind for breach of warranty or negligence.