

LACHAT INSTRUMENT'S OMNION 4.0 SOFTWARE

A REVOLUTION IN DATA ANALYSIS SOFTWARE!

INTRODUCTION

Omnion Software from Lachat Instruments is the first software to offer a completely unified computing environment for **both** flow injection analysis and ion chromatography! Easy-to-use, but tremendously powerful and versatile, Omnion 4.0 simplifies the operation of the QuikChem® 8500 Series Automated Ion Analyzers.



WINDOWS 7, WINDOWS 8.X AND WINDOWS 10 COMPATIBLE

- This 64-bit software system is entirely compatible with Windows 7, Windows 8.X and Windows 10.
- Saves time and money by delivering a product that is compatible with the most current IT systems.

INTEGRATED RUN WORKSHEET

- User interface saves training time and reduces opportunities for error.
- Saves operator time by making intelligent run-time decisions, and monitoring run progress.
- Integrated Data Quality Management (DQM) to Schedule, Identify, and Report Quality Assurance Samples.
- DQM functions include: Check Standards, Duplicates, Spike, Spiked Duplicates, Calibration Standard, Blank, Lab Control Sample
- The Omnion software interface is provided in multiple languages - including Spanish, German, French, Portuguese, and Italian

SEAMLESS DATA HANDLING


- Creates data that seamlessly exports to LIMS.
- Eliminates or minimizes the need for custom data export/review programs.
- Saves time during review and transfer of data.
- Custom Reports

SAMPLE FLAGGING CAPABILITIES

- Samples with pre-defined criteria are easy to find, saving report reviewing time by operator and QA staff.
- Saves data analysis time by eliminating unneeded or unwanted data in exported sample reports.
- Saves personnel time by creating more concise reports for import into LIMS.

SIMPLE RE-RUNS

- Saves time by eliminating operator errors during startup and minimizes need for a manual or technical support.



Specification	Minimum	Recommended
• Processor	• 1GHz	• Dual Core 1.8GHz (or better)
• Memory	• 2GB RAM	• 4GB RAM (or better)
• Hard Disk Free Space	• 500 MB	• 3000 MB (or better)
• Display	• 1024 x 768 resolution	• 1080 x 1024 resolution
• 17" Flat Panel Monitor	• 17" Flat Panel Monitor	
• Operating System	• 64-bit Windows 7, Windows 8.X, or Windows 10	
• Interfaces	• CD/DVD-ROM, USB 1.1 or 2.0 (x2) for instrument and autosampler connections	

How To Order OMNION 4.0 SOFTWARE:

PN: **A85081** Omnion 4.0 Software CD - FIA/IC

In the U.S.A.
 Call toll-free (800) 247-7613.

Outside the U.S.A.
 Contact the LACHAT distributor serving you.
www.lachatinstruments.com or www.hach.com

Email:
sales@lachatinstruments.com or techhelp@hach.com

LL011 Rev. 2

Lachat Instruments Brand Loveland, Colorado USA

United States: 800-247-7613 tel | 970-669-3050 tel | 970-461-3915 fax | sales@lachatinstruments.com | techhelp@hach.com
 Outside United States: Contact the Lachat office or distributor serving you. | sales@lachatinstruments.com | intltech@hach.com
www.lachatinstruments.com

Printed in U.S.A. ©Hach Company, 2016. All rights reserved.

